



CONTROL DM-RET

DM-RET

HEMATOLOGY CONTROL

CONTROL


 2022-12-11

 2023-03-10


Applicable Instruments	Parameter	Unit	LOT	RET0123L	LOT	RET0123N	LOT	RET0123H
DYMIND DH-600 DH-602 DH-605 DH-610 DH-612 DH-615 (Technical File Version A1.0 or higher)	RET#	$\times 10^{12}/L$	0.0085	± 0.0080	0.2634	± 0.1316	0.3414	± 0.1200
	RET%	%	0.15	± 0.15	5.48	± 2.00	13.45	± 4.50
	LFR	%	89.5	± 27.0	90.0	± 27.0	90.0	± 27.0
	MFR	%	8.0	± 12.0	8.0	± 12.0	8.0	± 12.0
	HFR	%	2.6	± 3.5	2.0	± 3.0	2.0	± 3.0
	IRF	%	10.5	± 15.0	10.1	± 15.0	10.0	± 15.0
	RHE	pg	14.9	± 8.0	16.7	± 6.0	15.2	± 6.0

【NOTE】


1. The controls should be stored in refrigerator (2°C~8°C). After opening, it will keep stable for 14 days when it is stored airtight at 2°C~8°C.
2. Before mixing and running the control after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15°C~30°C).
3. Controls must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
4. After using, put the controls back into the refrigerator to prevent contamination and evaporation.



RET0123L



RET0123N



RET0123H