	文件类型	技术文件	文件版本	V3.4.0
	项目代号	/	页 码	第 0 页 共 1 页
	文件编号	GL DM-56-18-en		

文件名称:

CONTROL DM-6D 6D0823 2023-10-10

适用范围: 1116 项目

编制: 夏琳

审核: 刘金盈

批准: 郑剑通

日期: 2023.08.02

日期: 2023.08.02

日期: 2023.08.02

修 订 记 录

版本	TCN/ECR/PCN 编号	修订内容概要	修改人	批准人	批准日期
V3.4.0	/	新建	夏琳	郑剑通	2023.08.02

发至:

总经办 质量部 财务部 营销部 客服 计划 采购 工程 仓库


设备研发 试剂研发 物料质量 设备生产 试剂生产 试剂质检 设备质检


由试剂研发

DM-6D

HEMATOLOGY CONTROL

Reference Values provided by DYMIND


CONTROL
 2023-07-13

 2023-10-10


Instruments	Parameter	Unit	LOT	6D0823L	LOT	6D0823N	LOT	6D0823H
DYMIND DH-600 DH-602 DH-605 DH-610 DH-612 DH-615 (Technical File Version A2.0 or higher)	WBC	×10 ⁹ /L	4.12	±1.20	7.69	±1.80	22.73	±3.00
	Neu#	×10 ⁹ /L	2.91	±0.40	4.66	±0.80	15.55	±2.00
	Lym#	×10 ⁹ /L	0.62	±0.40	2.16	±0.70	5.08	±1.80
	Mon#	×10 ⁹ /L	0.34	±0.40	0.49	±0.50	0.88	±0.80
	Eos#	×10 ⁹ /L	0.25	±0.40	0.38	±0.70	1.22	±1.60
	Bas#	×10 ⁹ /L	3.29	±0.40	5.46	±0.70	16.48	±1.60
	IG#	×10 ⁹ /L	0.52	±0.50	0.82	±1.00	2.65	±2.50
	Neu%	%	70.6	±12.0	60.6	±12.0	68.4	±12.0
	Lym%	%	15.0	±10.0	28.1	±10.0	22.4	±10.0
	Mon%	%	8.4	±8.0	6.3	±8.0	3.9	±8.0
	Eos%	%	6.1	±10.0	4.9	±10.0	5.4	±10.0
	Bas%	%	79.8	±10.0	71.0	±10.0	72.5	±10.0
	IG%	%	12.5	±10.0	10.7	±10.0	11.6	±10.0
	RBC	×10 ¹² /L	2.33	±0.20	4.45	±0.30	5.52	±0.50
	HGB	g/L	61	±4	134	±6	176	±8
	HCT	%	18.1	±2.0	39.5	±2.5	51.3	±3.0
	MCV	fL	77.7	±5.0	88.8	±5.0	92.9	±5.0
	MCH	pg	26.4	±2.5	30.1	±2.5	31.9	±2.5
	MCHC	g/L	339	±30	339	±30	343	±30
	RDW-CV	%	16.5	±5.0	14.6	±5.0	14.1	±6.0
	RDW-SD	fL	46.9	±10.0	47.9	±10.0	48.0	±12.0
	PLT	×10 ⁹ /L	76	±20	289	±45	553	±65
	MPV	fL	7.5	±3.0	7.7	±3.0	7.5	±3.0
	PDW	fL	9.5	±5.0	9.3	±5.0	9.7	±5.0
	PCT	%	0.057	±0.050	0.221	±0.100	0.417	±0.200
	P-LCR	%	12.4	±20.0	12.3	±20.0	11.5	±20.0
	P-LCC	×10 ⁹ /L	10	±10	35	±30	64	±50
IPF	%	0.8	±0.8	0.5	±0.5	0.5	±0.5	

【NOTE】


1. The controls should be stored in refrigerator (2°C~8°C). After opening, it will keep stable for 14 days when it is stored airtight at 2°C~8°C.
2. Before mixing and running the control after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15°C~30°C).
3. Controls must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
4. After using, put the controls back into the refrigerator to prevent contamination and evaporation.



6D0823L



6D0823N



6D0823H