

## Table of CRP Calibrator Reference Value

### C-Reactive Protein Calibrator

**Code: DM3113**

**LOT** 2022122001

 20231219

Model	Level (multi-point)	Reference Value (mg/L)
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.9)	<b>Calibrator 1</b>	5.4
	<b>Calibrator 2</b>	20.5
	<b>Calibrator 3</b>	48.6
	<b>Calibrator 4</b>	88.9
	<b>Calibrator 5</b>	134.9
	<b>Calibrator 6</b>	175.0
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above with heating)	<b>Calibrator 1</b>	4.0
	<b>Calibrator 2</b>	19.1
	<b>Calibrator 3</b>	50.3
	<b>Calibrator 4</b>	105.6
	<b>Calibrator 5</b>	147.2
	<b>Calibrator 6</b>	186.6
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above without heating)	<b>Calibrator 1</b>	5.3
	<b>Calibrator 2</b>	22.1
	<b>Calibrator 3</b>	50.7
	<b>Calibrator 4</b>	104.0
	<b>Calibrator 5</b>	152.6
	<b>Calibrator 6</b>	213.7
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A5.11 and below)	<b>Calibrator 1</b>	4.1
	<b>Calibrator 2</b>	20.7
	<b>Calibrator 3</b>	50.1
	<b>Calibrator 4</b>	99.3
	<b>Calibrator 5</b>	146.8
	<b>Calibrator 6</b>	184.4
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above with heating)	<b>Calibrator 1</b>	4.8
	<b>Calibrator 2</b>	22.7
	<b>Calibrator 3</b>	51.3
	<b>Calibrator 4</b>	103.2
	<b>Calibrator 5</b>	154.3
	<b>Calibrator 6</b>	194.0

DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above without heating)	<b>Calibrator 1</b>	5.5
	<b>Calibrator 2</b>	23.7
	<b>Calibrator 3</b>	53.5
	<b>Calibrator 4</b>	109.6
	<b>Calibrator 5</b>	158.5
	<b>Calibrator 6</b>	207.3
D6-CRP, DH76CRP (Applicable to technical document version B2.0)	<b>Calibrator 1</b>	4.7
	<b>Calibrator 2</b>	21.1
	<b>Calibrator 3</b>	54.5
	<b>Calibrator 4</b>	100.7
	<b>Calibrator 5</b>	156.2
	<b>Calibrator 6</b>	200.4
D6-CRP, DH76CRP (Applicable to technical document version B2.2/B2.3/B2.4/B2.5)	<b>Calibrator 1</b>	4.0
	<b>Calibrator 2</b>	19.5
	<b>Calibrator 3</b>	51.5
	<b>Calibrator 4</b>	105.3
	<b>Calibrator 5</b>	147.0
	<b>Calibrator 6</b>	190.8
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.4 and below with heating)	<b>Calibrator 1</b>	4.8
	<b>Calibrator 2</b>	20.1
	<b>Calibrator 3</b>	53.4
	<b>Calibrator 4</b>	105.4
	<b>Calibrator 5</b>	156.2
	<b>Calibrator 6</b>	205.0
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.4 and below without heating)	<b>Calibrator 1</b>	5.3
	<b>Calibrator 2</b>	22.3
	<b>Calibrator 3</b>	54.7
	<b>Calibrator 4</b>	105.2
	<b>Calibrator 5</b>	154.8
	<b>Calibrator 6</b>	200.5
DM79X, DM78X, DM72X (Applicable to technical document version A1.0 and above ) DH36X, DH33X, DH31X (Applicable to technical document versions A1.0 and above)	<b>Calibrator 1</b>	4.9
	<b>Calibrator 2</b>	20.5
	<b>Calibrator 3</b>	55.4
	<b>Calibrator 4</b>	106.5
	<b>Calibrator 5</b>	156.0
	<b>Calibrator 6</b>	203.1

DM79X、DM78X、DM72X (Applicable to technical document version A1.0 and above )	<b>Calibrator 1</b>	4.3
	<b>Calibrator 2</b>	21.5
	<b>Calibrator 3</b>	57.0
	<b>Calibrator 4</b>	107.8
	<b>Calibrator 5</b>	147.2
	<b>Calibrator 6</b>	180.3