



Table of CRP Control Reference Value

C-Reactive Protein Control

Code: DM3513

 2023011001

 20240109

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DH-100T、DH-102T、DH-103T、 DH-105T、DH-106T、DH-108T、 DH-2000T[H1T1]、 DH-2000T[H2T1]、 DH-2000T[H1T2]、 DH-2000T[H2T2] (Applicable to technical document version A1.0 and above)	Control 1	2.7	1.9-3.5	±0.8
	Control 2	17.3	12.1-22.5	±5.2
	Control 3	50.7	35.5-65.9	±15.2