


Table of CRP Calibrator Reference Value

C-Reactive Protein Calibrator

Code: DM3513

LOT 2023011001

 20240109

Model	Level	Reference Value (mg/L)
DH-100T、DH-102T、DH-103T、DH-105T、DH-106T、 DH-108T、DH-2000T[H1T1]、DH-2000T[H2T1]、 DH-2000T[H1T2]、DH-2000T[H2T2] (Applicable to technical document version A1.0 and above)	Calibrator 1	5.5
	Calibrator 2	25.3
	Calibrator 3	64.2
	Calibrator 4	121.6
	Calibrator 5	168.2
	Calibrator 6	235.1