

	文件类型	技术文件	文件版本	V3.5.2	
	项目代号	DM3113	页 码	第 0 页 共 7 页	
	文件编号	GL DM-56-08-en			
文件名称: <p style="text-align: center;">Table of CRP Control Reference Value</p>					
适用范围: <p style="text-align: center;">DM3113 项目</p>					
编制: 曾婷婷 日期: 2023.02.27		审核: 曾婷婷 日期: 2023.02.27		批准: 秦世界 日期: 2023.02.27	
修 订 记 录					
版本	TCN/ECR/P CN 编号	修订内容概要	修改人	批准人	批准日期
V3.5.0	/	新建	秦俊飞	秦世界	2022.12.30
V3.5.1	/	更新 1109 加热模式靶值	曾婷婷	秦世界	2023.01.11
V3.5.2	/	同机型分页, 更新 P101 靶值	曾婷婷	魏保立	2023.02.27
发至:	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input type="checkbox"/> 营销部 <input type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input type="checkbox"/> 设备研发 <input type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input type="checkbox"/> 试剂质检 <input type="checkbox"/> 设备质 检				
由试剂研发					

Table of CRP Control Reference Value

C-Reactive Protein Control

Code: DM3113

LOT 2022122001

🕒 20231219

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.9)	Control 1	2.7	1.9-3.5	±0.8
	Control 2	15.0	10.5-19.5	±4.5
	Control 3	40.5	28.3-52.7	±12.2
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above with heating)	Control 1	2.6	1.8-3.4	±0.8
	Control 2	13.8	9.7-17.9	±4.1
	Control 3	40.3	28.2-52.4	±12.1
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above without heating)	Control 1	2.9	2.0-3.8	±0.9
	Control 2	15.6	10.9-20.3	±4.7
	Control 3	41.3	28.9-53.7	±12.4
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A5.11 and below)	Control 1	2.3	1.6-3.0	±0.7
	Control 2	14.4	10.1-18.7	±4.3
	Control 3	38.6	27.0-50.2	±11.6
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above with heating)	Control 1	2.8	2.0-3.6	±0.8
	Control 2	15.7	11.0-20.4	±4.7
	Control 3	42.6	29.8-55.4	±12.8
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above without heating)	Control 1	3.1	2.2-4.0	±0.9
	Control 2	16.5	11.5-21.5	±5.0
	Control 3	44.6	31.2-58.0	±13.4

D6-CRP, DH76CRP (Applicable to technical document version B2.0)	Control 1	3.0	2.1-3.9	± 0.9
	Control 2	14.9	10.4-19.4	± 4.5
	Control 3	37.3	26.1-48.5	± 11.2
D6-CRP, DH76CRP (Applicable to technical document version B2.2/B2.3/B2.4/B2.5)	Control 1	2.2	1.5-2.9	± 0.7
	Control 2	13.9	9.7-18.1	± 4.2
	Control 3	42.7	29.9-55.5	± 12.8
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.4 and below with heating)	Control 1	2.7	1.9-3.5	± 0.8
	Control 2	16.0	11.2-20.8	± 4.8
	Control 3	39.7	27.8-51.6	± 11.9
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.4 and below without heating)	Control 1	3.0	2.1-3.9	± 0.9
	Control 2	15.2	10.6-19.8	± 4.6
	Control 3	41.3	28.9-53.7	± 12.4
DM79X, DM78X, DM72X (Applicable to technical document version A1.0 and above)	Control 1	3.1	2.2-4.0	± 0.9
	Control 2	14.9	10.4-19.4	± 4.5
	Control 3	42.4	29.7-55.1	± 12.7
DH36X, DH33X, DH31X (Applicable to technical document versions A1.0 and above)	Control 1	2.5	1.7-3.3	± 0.8
	Control 2	13.4	9.4-17.4	± 4.0
	Control 3	43.5	30.4-56.6	± 13.1
DP-H10, DP-H12, DP-H13, DP-H15, DP-H16, DP-H17 (Applicable to technical document versions A1.7 and above)	Control 1	3.3	2.3-4.3	± 1.0
	Control 2	16.5	11.5-21.5	± 5.0
	Control 3	44.1	30.9-57.3	± 13.2