
	文件类型	技术文件	文件版本	V3.2.0		
	项目代号	/	页 码	第0页 共 1页		
	文件编号	GL DM-56-17-en				
文件名称： <p style="text-align: center;">CONTROL DM-RET RET0323 2023-05-10</p> 适用范围：1116 项目						
编制：		审核：		批准：		
日期：		日期：		日期：		
修 订 记 录						
版本	TCN/ECR/PCN 编号	修订内容概要		修改人	批准人	批准日期
V3.2.0	/	新建		段利超	郑剑通	2023.03.13
发至：	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input checked="" type="checkbox"/> 营销部 <input checked="" type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input checked="" type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input checked="" type="checkbox"/> 设备研发 <input checked="" type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input checked="" type="checkbox"/> 试剂质检 <input checked="" type="checkbox"/> 设备质检					
由试剂研发						

CONTROL DM-RET

DM-RET

HEMATOLOGY CONTROL


CONTROL
 2023-02-10

 2023-05-10


Applicable Instruments	Parameter	Unit	LOT RET0323L	LOT RET0323N	LOT RET0323H
DYMIND DH-600 DH-602 DH-605 DH-610 DH-612 DH-615 (Technical File Version A1.6 or higher)	RET#	$\times 10^{12}/L$	0.0199 ±0.0100	0.3014 ±0.1316	0.3352 ±0.1200
	RET%	%	0.35 ±0.15	6.49 ±2.00	14.19 ±4.50
	LFR	%	89.8 ±27.0	90.0 ±27.0	89.9 ±27.0
	MFR	%	8.0 ±12.0	8.0 ±12.0	8.1 ±12.0
	HFR	%	2.2 ±3.5	2.0 ±3.0	2.0 ±3.0
	IRF	%	10.2 ±15.0	10.0 ±15.0	10.1 ±15.0
	RHE	pg	15.6 ±8.0	14.3 ±6.0	15.8 ±8.0

【NOTE】


1. The controls should be stored in refrigerator (2°C~8°C). After opening, it will keep stable for 14 days when it is stored airtight at 2°C~8°C.
2. Before mixing and running the control after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15°C~30°C).
3. Controls must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
4. After using, put the controls back into the refrigerator to prevent contamination and evaporation.



RET0323L



RET0323N



RET0323H