

	文件类型	技术文件	文件版本	V4.1.0	
	项目代号	DM3113	页 码	第 0 页 共 2 页	
	文件编号	GL DM-56-08-en			
文件名称： <p style="text-align: center;">Table of CRP Control Reference Value</p>					
适用范围： <p style="text-align: center;">DM3113 项目</p>					
编制： 曾婷婷 日期： 2023.03.28		审核： 魏保立 日期： 2023.03.28		批准： 魏保立 日期： 2023.03.28	
修 订 记 录					
版本	TCN/ECR/P CN 编号	修订内容概要	修改人	批准人	批准日期
V4.1.0	/	新建	曾婷婷	魏保立	2023.03.28
发至：	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input type="checkbox"/> 营销部 <input type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input type="checkbox"/> 设备研发 <input type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input type="checkbox"/> 试剂质检 <input type="checkbox"/> 设备质 检				
由试剂研发					

Table of CRP Control Reference Value

C-Reactive Protein Control

Code: DM3113

LOT 2023030201

🕒 20240301

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.9)	Control 1	2.2	1.5-2.9	±0.7
	Control 2	16.2	11.3-21.1	±4.9
	Control 3	34.9	24.4-45.4	±10.5
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above with heating)	Control 1	2.7	1.9-3.5	±0.8
	Control 2	14.3	10.0-18.6	±4.3
	Control 3	33.4	23.4-43.4	±10.0
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above without heating)	Control 1	3.2	2.2-4.2	±1.0
	Control 2	15.3	10.7-19.9	±4.6
	Control 3	37.1	26.0-48.2	±11.1
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A5.11 and below)	Control 1	2.0	1.4-2.6	±0.6
	Control 2	14.7	10.3-19.1	±4.4
	Control 3	32.9	23.0-42.8	±9.9
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above with heating)	Control 1	2.6	1.8-3.4	±0.8
	Control 2	15.2	10.6-19.8	±4.6
	Control 3	37.7	26.4-49.0	±11.3
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above without heating)	Control 1	2.9	2.0-3.8	±0.9
	Control 2	16.1	11.3-20.9	±4.8
	Control 3	37.4	26.2-48.6	±11.2

D6-CRP, DH76CRP (Applicable to technical document version B2.0)	Control 1	2.6	1.8-3.4	± 0.8
	Control 2	14.5	10.1-18.9	± 4.4
	Control 3	34.9	24.4-45.4	± 10.5
D6-CRP, DH76CRP (Applicable to technical document version B2.2 and above)	Control 1	2.9	2.0-3.8	± 0.9
	Control 2	15.7	11.0-20.4	± 4.7
	Control 3	37.6	26.3-48.9	± 11.3
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.5 and below with heating)	Control 1	3.1	2.2-4.0	± 0.9
	Control 2	15.9	11.1-20.7	± 4.8
	Control 3	32.4	22.7-42.1	± 9.7
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.5 and below without heating)	Control 1	3.2	2.2-4.2	± 1.0
	Control 2	16.0	11.2-20.8	± 4.8
	Control 3	38.2	26.7-49.7	± 11.5
DM79X, DM78X, DM72X (Applicable to technical document version A1.0 and above)	Control 1	3.1	2.2-4.0	± 0.9
	Control 2	16.6	11.6-21.6	± 5.0
	Control 3	37.3	26.1-48.5	± 11.2
DH36X, DH33X, DH31X (Applicable to technical document versions A1.0 and above)	Control 1	3.2	2.2-4.2	± 1.0
	Control 2	15.3	10.7-19.9	± 4.6
	Control 3	34.9	24.4-45.4	± 10.5
DP-H10, DP-H12, DP-H13, DP-H15, DP-H16, DP-H17 (Applicable to technical document versions A1.7 and above)	Control 1	3.5	2.5-4.6	± 1.1
	Control 2	18.2	12.7-23.7	± 5.5
	Control 3	39.7	27.8-51.6	± 11.9