

	文件类型	技术文件	文件版本	V4.1.0	
	项目代号	DM3113	页 码	第 0 页 共 3 页	
	文件编号	GL DM-56-09-en			
文件名称: <p style="text-align: center;"><b>Table of CRP Calibrator Reference Value</b></p>					
适用范围: <p style="text-align: center;">DM3113 项目</p>					
编制: 曾婷婷 日期 2023.03.28		审核: 魏保立 日期: 2023.03.28		批准: 魏保立 日期: 2023.03.28	
修 订 记 录					
版本	TCN/ECR/P CN 编号	修订内容概要	修改人	批准人	批准日期
V4.1.0	/	新建	曾婷婷	魏保立	2023.03.28
发至:	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input type="checkbox"/> 营销部 <input type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input type="checkbox"/> 设备研发 <input type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input type="checkbox"/> 试剂质检 <input type="checkbox"/> 设备质 检				
由试剂研发					

## Table of CRP Calibrator Reference Value

### C-Reactive Protein Calibrator

**Code: DM3113**

LOT 2023030201

 20240301

Model	Level (multi-point)	Reference Value (mg/L)
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.9)	Calibrator 1	4.8
	Calibrator 2	18.7
	Calibrator 3	44.1
	Calibrator 4	88.9
	Calibrator 5	134.9
	Calibrator 6	175.0
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above with heating)	Calibrator 1	4.8
	Calibrator 2	19.1
	Calibrator 3	44.1
	Calibrator 4	95.3
	Calibrator 5	146.8
	Calibrator 6	190.5
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above without heating)	Calibrator 1	5.5
	Calibrator 2	20.0
	Calibrator 3	45.7
	Calibrator 4	102.6
	Calibrator 5	145.7
	Calibrator 6	190.3
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A5.11 and below)	Calibrator 1	4.5
	Calibrator 2	18.8
	Calibrator 3	46.3
	Calibrator 4	94.1
	Calibrator 5	127.8
	Calibrator 6	150.3
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above with heating)	Calibrator 1	4.8
	Calibrator 2	19.4
	Calibrator 3	51.3
	Calibrator 4	103.2
	Calibrator 5	154.3
	Calibrator 6	194.6

DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above without heating)	<b>Calibrator 1</b>	5.4
	<b>Calibrator 2</b>	20.9
	<b>Calibrator 3</b>	46.4
	<b>Calibrator 4</b>	98.4
	<b>Calibrator 5</b>	142.5
	<b>Calibrator 6</b>	188.9
D6-CRP, DH76CRP (Applicable to technical document version B2.0)	<b>Calibrator 1</b>	4.5
	<b>Calibrator 2</b>	19.2
	<b>Calibrator 3</b>	43.6
	<b>Calibrator 4</b>	100.4
	<b>Calibrator 5</b>	141.9
	<b>Calibrator 6</b>	200.4
D6-CRP, DH76CRP (Applicable to technical document version B2.2 and above)	<b>Calibrator 1</b>	5.2
	<b>Calibrator 2</b>	20.3
	<b>Calibrator 3</b>	46.5
	<b>Calibrator 4</b>	101.5
	<b>Calibrator 5</b>	143.7
	<b>Calibrator 6</b>	172.2
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.5 and below with heating)	<b>Calibrator 1</b>	4.9
	<b>Calibrator 2</b>	20.3
	<b>Calibrator 3</b>	45.3
	<b>Calibrator 4</b>	98.5
	<b>Calibrator 5</b>	153.4
	<b>Calibrator 6</b>	188.7
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.5 and below without heating)	<b>Calibrator 1</b>	5.3
	<b>Calibrator 2</b>	20.4
	<b>Calibrator 3</b>	48.6
	<b>Calibrator 4</b>	98.9
	<b>Calibrator 5</b>	149.7
	<b>Calibrator 6</b>	194.5
DM79X, DM78X, DM72X (Applicable to technical document version A1.0 and above)	<b>Calibrator 1</b>	5.4
	<b>Calibrator 2</b>	20.6
	<b>Calibrator 3</b>	45.4
	<b>Calibrator 4</b>	99.1
	<b>Calibrator 5</b>	139.5
	<b>Calibrator 6</b>	171.2

DH36X、DH33X、DH31X (Applicable to technical document versions A1.0 and above)	<b>Calibrator 1</b>	4.9
	<b>Calibrator 2</b>	19.8
	<b>Calibrator 3</b>	45.6
	<b>Calibrator 4</b>	99.7
	<b>Calibrator 5</b>	145.3
	<b>Calibrator 6</b>	174.6