


Table of CRP Control Reference Value**C-Reactive Protein Control****Code: DM3513**

LOT 2023021701

 20240216

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DH-100T、DH-102T、DH-103T、 DH-105T、DH-106T、DH-108T、 DH-2000T[H1T1]、 DH-2000T[H2T1]、 DH-2000T[H1T2]、 DH-2000T[H2T2] (Applicable to technical document version A1.0 and above)	Control 1	2.8	2.0-3.6	±0.8
	Control 2	19.4	13.6-25.2	±5.8
	Control 3	51.1	35.8-66.4	±15.3