



## Table of CRP Control Reference Value

### C-Reactive Protein Control

**Code: DM3113**

LOT 2023030201

🕒 20240301

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.9)	<b>Control 1</b>	2.2	1.5-2.9	±0.7
	<b>Control 2</b>	16.2	11.3-21.1	±4.9
	<b>Control 3</b>	34.9	24.4-45.4	±10.5
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above with heating)	<b>Control 1</b>	2.7	1.9-3.5	±0.8
	<b>Control 2</b>	14.3	10.0-18.6	±4.3
	<b>Control 3</b>	33.4	23.4-43.4	±10.0
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above without heating)	<b>Control 1</b>	3.2	2.2-4.2	±1.0
	<b>Control 2</b>	15.3	10.7-19.9	±4.6
	<b>Control 3</b>	37.1	26.0-48.2	±11.1
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A5.11 and below)	<b>Control 1</b>	2.0	1.4-2.6	±0.6
	<b>Control 2</b>	14.7	10.3-19.1	±4.4
	<b>Control 3</b>	32.9	23.0-42.8	±9.9
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above with heating)	<b>Control 1</b>	2.6	1.8-3.4	±0.8
	<b>Control 2</b>	15.2	10.6-19.8	±4.6
	<b>Control 3</b>	37.7	26.4-49.0	±11.3
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above without heating)	<b>Control 1</b>	2.9	2.0-3.8	±0.9
	<b>Control 2</b>	16.1	11.3-20.9	±4.8
	<b>Control 3</b>	37.4	26.2-48.6	±11.2

D6-CRP, DH76CRP (Applicable to technical document version B2.0)	<b>Control 1</b>	2.6	1.8-3.4	$\pm 0.8$
	<b>Control 2</b>	14.5	10.1-18.9	$\pm 4.4$
	<b>Control 3</b>	34.9	24.4-45.4	$\pm 10.5$
D6-CRP, DH76CRP (Applicable to technical document version B2.2 and above)	<b>Control 1</b>	2.9	2.0-3.8	$\pm 0.9$
	<b>Control 2</b>	15.7	11.0-20.4	$\pm 4.7$
	<b>Control 3</b>	37.6	26.3-48.9	$\pm 11.3$
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.5 and below with heating)	<b>Control 1</b>	3.1	2.2-4.0	$\pm 0.9$
	<b>Control 2</b>	15.9	11.1-20.7	$\pm 4.8$
	<b>Control 3</b>	35.7	25.0-46.4	$\pm 10.7$
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.5 and below without heating)	<b>Control 1</b>	3.2	2.2-4.2	$\pm 1.0$
	<b>Control 2</b>	16.0	11.2-20.8	$\pm 4.8$
	<b>Control 3</b>	38.2	26.7-49.7	$\pm 11.5$
DM79X, DM78X, DM72X (Applicable to technical document version A1.0 and above)	<b>Control 1</b>	3.1	2.2-4.0	$\pm 0.9$
	<b>Control 2</b>	16.6	11.6-21.6	$\pm 5.0$
	<b>Control 3</b>	37.3	26.1-48.5	$\pm 11.2$
DH36X, DH33X, DH31X (Applicable to technical document versions A1.0 and above)	<b>Control 1</b>	3.2	2.2-4.2	$\pm 1.0$
	<b>Control 2</b>	15.3	10.7-19.9	$\pm 4.6$
	<b>Control 3</b>	34.9	24.4-45.4	$\pm 10.5$
DP-H10, DP-H12, DP-H13, DP-H15, DP-H16, DP-H17 (Applicable to technical document versions A1.7 and above)	<b>Control 1</b>	3.5	2.5-4.6	$\pm 1.1$
	<b>Control 2</b>	18.2	12.7-23.7	$\pm 5.5$
	<b>Control 3</b>	39.7	27.8-51.6	$\pm 11.9$