

	文件类型	技术文件	文件版本	V4.3.0	
	项目代号	DM3113	页 码	第 0 页 共 2 页	
	文件编号	GL DM-56-08-en			
文件名称: <p style="text-align: center;"><b>Table of CRP Control Reference Value</b></p>					
适用范围: <p style="text-align: center;">DM3113 项目</p>					
编制: 曾婷婷 日期: 2023.06.08		审核: 魏保立 日期: 2023.06.08		批准: 魏保立 日期: 2023.06.08	
修 订 记 录					
版本	TCN/ECR/P CN 编号	修订内容概要	修改人	批准人	批准日期
V4.3.0	/	新建	曾婷婷	魏保立	2023.06.08
发至:	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input type="checkbox"/> 营销部 <input type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input type="checkbox"/> 设备研发 <input type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input type="checkbox"/> 试剂质检 <input type="checkbox"/> 设备质 检				
由试剂研发					

## Table of CRP Control Reference Value

### C-Reactive Protein Control

**Code: DM3113**

LOT 2023052601

🕒 20240525

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.9)	<b>Control 1</b>	2.8	2.0-3.6	±0.8
	<b>Control 2</b>	16.5	11.5-21.5	±5.0
	<b>Control 3</b>	37.7	26.4-49.0	±11.3
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above with heating)	<b>Control 1</b>	4.0	2.8-5.2	±1.2
	<b>Control 2</b>	17.3	12.1-22.5	±5.2
	<b>Control 3</b>	39.6	27.7-51.5	±11.9
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above without heating)	<b>Control 1</b>	3.9	2.7-5.1	±1.2
	<b>Control 2</b>	17.8	12.5-23.1	±5.3
	<b>Control 3</b>	43.8	30.7-56.9	±13.1
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A5.11 and below)	<b>Control 1</b>	2.9	2.0-3.8	±0.9
	<b>Control 2</b>	15.5	10.8-20.2	±4.7
	<b>Control 3</b>	40.0	28.0-52.0	±12.0
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above with heating)	<b>Control 1</b>	3.3	2.3-4.3	±1.0
	<b>Control 2</b>	17.4	12.2-22.6	±5.2
	<b>Control 3</b>	44.2	30.9-57.5	±13.3
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above without heating)	<b>Control 1</b>	4.4	3.1-5.7	±1.3
	<b>Control 2</b>	19.0	13.3-24.7	±5.7
	<b>Control 3</b>	43.2	30.2-56.2	±13.0

D6-CRP, DH76CRP (Applicable to technical document version B2.0)	<b>Control 1</b>	3.8	2.7-4.9	$\pm 1.1$
	<b>Control 2</b>	17.7	12.4-23.0	$\pm 5.3$
	<b>Control 3</b>	41.5	29.0-54.0	$\pm 12.5$
D6-CRP, DH76CRP (Applicable to technical document version B2.2 and above)	<b>Control 1</b>	3.3	2.3-4.3	$\pm 1.0$
	<b>Control 2</b>	16.7	11.7-21.7	$\pm 5.0$
	<b>Control 3</b>	40.2	28.1-52.3	$\pm 12.1$
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.5 and below with heating)	<b>Control 1</b>	4.0	2.8-5.2	$\pm 1.2$
	<b>Control 2</b>	18.5	12.9-24.1	$\pm 5.6$
	<b>Control 3</b>	43.5	30.4-56.6	$\pm 13.1$
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.5 and below without heating)	<b>Control 1</b>	5.0	3.5-6.5	$\pm 1.5$
	<b>Control 2</b>	19.9	13.9-25.9	$\pm 6.0$
	<b>Control 3</b>	43.4	30.4-56.4	$\pm 13.0$
DM79X, DM78X, DM72X (Applicable to technical document version A1.0 and above)	<b>Control 1</b>	4.1	2.9-5.3	$\pm 1.2$
	<b>Control 2</b>	18.0	12.6-23.4	$\pm 5.4$
	<b>Control 3</b>	41.4	29.0-53.8	$\pm 12.4$
DH36X, DH33X, DH31X (Applicable to technical document versions A1.0 and above)	<b>Control 1</b>	4.2	2.9-5.5	$\pm 1.3$
	<b>Control 2</b>	18.4	12.9-23.9	$\pm 5.5$
	<b>Control 3</b>	42.1	29.5-54.7	$\pm 12.6$
DP-H10, DP-H12, DP-H13, DP-H15, DP-H16, DP-H17 (Applicable to technical document versions A1.7 and above)	<b>Control 1</b>	/	/	/
	<b>Control 2</b>	/	/	/
	<b>Control 3</b>	/	/	/