





# DM-CAL

## HEMATOLOGY CALIBRATOR

LOT N0923

 2023-08-12 2023-10-10

Applicable Instruments	Parameter	Unit	Reference Value
<b>DYMIND</b> <b>DP-H10、</b> <b>DP-H12、</b> <b>DP-H13、</b> <b>DP-H15、</b> <b>DP-H16、</b> <b>DP-H17</b>  (Technical File Version A1.10 or higher) ;	<b>WBC</b>	×10 <sup>9</sup> /L	9.74
	<b>RBC</b>	×10 <sup>12</sup> /L	4.88
	<b>HGB</b>	g/L	137
	<b>MCV</b>	fL	85.3
	<b>PLT</b>	×10 <sup>9</sup> /L	256

**【NOTE】**

1. The calibrator should be stored in refrigerator (2℃~8℃). After opening, it will keep stable for 7 days when it is stored airtight at 2℃~8℃.
2. Before mixing and running the calibrator after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15℃~30℃).
3. Calibrator must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
4. After using, put the calibrator back into the refrigerator to prevent contamination and evaporation.