




**Table of SAA Calibrator Reference Value****Serum Amyloid A Protein Calibrator**

LOT 2023081101

 20240810

<b>Model</b>	<b>Level (multi-point)</b>	<b>Reference Value (mg/L)</b>
DM79X、DM78X、DM72X (Applicable to technical document version A1.0 and above)	<b>Calibrator 1</b>	12.1
	<b>Calibrator 2</b>	21.7
	<b>Calibrator 3</b>	52.7
	<b>Calibrator 4</b>	106.7
	<b>Calibrator 5</b>	212.1
	<b>Calibrator 6</b>	260.4
	<b>Calibrator 7</b>	333.3
DP-H10、DP-H12、DP-H13、DP-H15、DP-H16、 DP-H17 (Applicable to technical document version A1.7 and above)	<b>Calibrator 1</b>	13.3
	<b>Calibrator 2</b>	23.6
	<b>Calibrator 3</b>	55.2
	<b>Calibrator 4</b>	112.3
	<b>Calibrator 5</b>	214.6