


	文件类型	技术文件	文件版本	V9.9.0		
	项目代号	/	页 码	第 0 页 共 1 页		
	文件编号	GL DM-56-02-en				
文件名称： <p style="text-align: center;">CALIBRATOR NK0923 2023-10-10 for DH56</p> 适用范围：1105 A4.1~A10.3、1101 A8.1 及以上						
编制：王雨婷 日期：2023.09.04		审核：刘金盈 日期：2023.09.04		批准：郑剑通 日期：2023.09.04		
修 订 记 录						
版本	TCN/ECR/PCN 编号	修订内容概要		修改人	批准人	批准日期
V9.9.0	/	新建		王雨婷	郑剑通	2023.09.04
发至：	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input checked="" type="checkbox"/> 营销部 <input checked="" type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input checked="" type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input checked="" type="checkbox"/> 设备研发 <input checked="" type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input checked="" type="checkbox"/> 试剂质检 <input checked="" type="checkbox"/> 设备质检					
由试剂研发						


DM-CAL PLUS

HEMATOLOGY CALIBRATOR

Reference Values provided by DYMIND
CALIBRATOR

LOT NK0923

 2023-08-12

 2023-10-10

Applicable Instruments	Parameter	Unit	Reference Value
DH71, DH76 (Technical File Version A4.1 to A10.3) DH51, DH53, DH56 (Technical File Version A8.1 or higher)	WBC	$\times 10^9/L$	9.57
	RBC	$\times 10^{12}/L$	4.54
	HGB	g/L	133
	MCV	fL	90.6
	PLT	$\times 10^9/L$	244

【NOTE】

1. The calibrator should be stored in refrigerator (2°C~8°C). After opening, it will keep stable for 7 days when it is stored airtight at 2°C~8°C.
2. Before mixing and running the calibrator after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15°C~30°C).
3. Calibrator must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
4. After using, put the calibrator back into the refrigerator to prevent contamination and evaporation.