

	文件类型	技术文件	文件版本	V4.5.0	
	项目代号	DM3113	页 码	第 0 页 共 2 页	
	文件编号	GL DM-56-08-en			
文件名称: <p style="text-align: center;">Table of CRP Control Reference Value</p>					
适用范围: <p style="text-align: center;">DM3113 项目</p>					
编制: 黄子苒 日期: 2023.10.16		审核: 曾婷婷 日期: 2023.10.16		批准: 魏保立 日期: 2023.10.16	
修 订 记 录					
版本	TCN/ECR/P CN 编号	修订内容概要	修改人	批准人	批准日期
V4.5.0	/	新建	黄子苒	魏保立	2023.10.16
发至:	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input type="checkbox"/> 营销部 <input type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input type="checkbox"/> 设备研发 <input type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input type="checkbox"/> 试剂质检 <input type="checkbox"/> 设备质 检				
由试剂研发					

Table of CRP Control Reference Value

C-Reactive Protein Control

Code: DM3113

LOT 2023092101

🕒 20240920

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.9)	Control 1	2.1	1.5-2.7	±0.6
	Control 2	11.4	8.0-14.8	±3.4
	Control 3	35.6	24.9-46.3	±10.7
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above with heating)	Control 1	2.8	2.0-3.6	±0.8
	Control 2	11.9	8.3-15.5	±3.6
	Control 3	35.9	25.1-46.7	±10.8
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above without heating)	Control 1	2.8	2.0-3.6	±0.8
	Control 2	12.3	8.6-16.0	±3.7
	Control 3	34.3	24.0-44.6	±10.3
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A5.11 and below)	Control 1	2.1	1.5-2.7	±0.6
	Control 2	11.7	8.2-15.2	±3.5
	Control 3	40.3	28.2-52.4	±12.1
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above with heating)	Control 1	3.0	2.1-3.9	±0.9
	Control 2	13.3	9.3-17.3	±4.0
	Control 3	37.4	26.2-48.6	±11.2
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above without heating)	Control 1	2.9	2.0-3.8	±0.9
	Control 2	13.2	9.2-17.2	±4.0
	Control 3	37.9	26.5-49.3	±11.4

D6-CRP, DH76CRP (Applicable to technical document version B2.0)	Control 1	2.7	1.9-3.5	± 0.8
	Control 2	12.5	8.8-16.3	± 3.8
	Control 3	38.7	27.1-50.3	± 11.6
D6-CRP, DH76CRP (Applicable to technical document version B2.2 and above)	Control 1	2.8	2.0-3.6	± 0.8
	Control 2	14.0	9.8-18.2	± 4.2
	Control 3	36.4	25.5-47.3	± 10.9
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.5 and below with heating)	Control 1	2.9	2.0-3.8	± 0.9
	Control 2	12.9	9.0-16.8	± 3.9
	Control 3	37.2	26.0-48.4	± 11.2
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.5 and below without heating)	Control 1	2.9	2.0-3.8	± 0.9
	Control 2	13.0	9.1-16.9	± 3.9
	Control 3	37.5	26.2-48.8	± 11.3
DM79X, DM78X, DM72X (Applicable to technical document version A1.0 and above)	Control 1	3.0	2.1-3.9	± 0.9
	Control 2	13.4	9.4-17.4	± 4.0
	Control 3	39.5	27.6-51.4	± 11.9
DH36X, DH33X, DH31X (Applicable to technical document versions A1.0 and above)	Control 1	3.0	2.1-3.9	± 0.9
	Control 2	12.7	8.9-16.5	± 3.8
	Control 3	38.1	26.7-49.5	± 11.4
DP-H10, DP-H12, DP-H13, DP-H15, DP-H16, DP-H17 (Applicable to technical document versions A1.7 and above)	Control 1	3.3	2.3-4.3	± 1.0
	Control 2	15.4	10.8-20.0	± 4.6
	Control 3	41.0	28.7-53.3	± 12.3