

Table of CRP Calibrator Reference Value

C-Reactive Protein Calibrator

Code: DM3113

LOT 2023092101

 20240920

Model	Level (multi-point)	Reference Value (mg/L)
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.9)	Calibrator 1	4.0
	Calibrator 2	16.4
	Calibrator 3	43.3
	Calibrator 4	79.9
	Calibrator 5	111.4
	Calibrator 6	139.6
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above with heating)	Calibrator 1	4.9
	Calibrator 2	18.1
	Calibrator 3	46.1
	Calibrator 4	93.2
	Calibrator 5	140.1
	Calibrator 6	185.4
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above without heating)	Calibrator 1	5.2
	Calibrator 2	17.9
	Calibrator 3	47.8
	Calibrator 4	102.6
	Calibrator 5	142.0
	Calibrator 6	185.3
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A5.11 and below)	Calibrator 1	4.2
	Calibrator 2	18.5
	Calibrator 3	52.6
	Calibrator 4	109.6
	Calibrator 5	154.3
	Calibrator 6	183.3
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above with heating)	Calibrator 1	5.2
	Calibrator 2	18.2
	Calibrator 3	46.9
	Calibrator 4	92.3
	Calibrator 5	139.8
	Calibrator 6	172.5

DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above without heating)	Calibrator 1	5.0
	Calibrator 2	17.9
	Calibrator 3	46.0
	Calibrator 4	94.8
	Calibrator 5	140.2
	Calibrator 6	172.4
D6-CRP, DH76CRP (Applicable to technical document version B2.0)	Calibrator 1	4.8
	Calibrator 2	18.5
	Calibrator 3	48.5
	Calibrator 4	93.5
	Calibrator 5	144.5
	Calibrator 6	187.9
D6-CRP, DH76CRP (Applicable to technical document version B2.2 and above)	Calibrator 1	5.2
	Calibrator 2	18.1
	Calibrator 3	46.2
	Calibrator 4	97.6
	Calibrator 5	135.3
	Calibrator 6	168.1
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.5 and below with heating)	Calibrator 1	5.2
	Calibrator 2	18.8
	Calibrator 3	48.1
	Calibrator 4	96.3
	Calibrator 5	145.5
	Calibrator 6	183.8
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A7.5 and below without heating)	Calibrator 1	5.1
	Calibrator 2	18.8
	Calibrator 3	47.1
	Calibrator 4	96.4
	Calibrator 5	147.7
	Calibrator 6	186.2
DM79X, DM78X, DM72X (Applicable to technical document version A1.0 and above)	Calibrator 1	5.3
	Calibrator 2	18.5
	Calibrator 3	50.7
	Calibrator 4	98.1
	Calibrator 5	143.3
	Calibrator 6	171.6

DH36X、DH33X、DH31X (Applicable to technical document versions A1.0 and above)	Calibrator 1	5.0
	Calibrator 2	18.3
	Calibrator 3	46.5
	Calibrator 4	92.5
	Calibrator 5	142.3
	Calibrator 6	181.6
DP-H10、DP-H12、DP-H13、DP-H15、DP-H16、 DP-H17 (Applicable to technical document versions A1.7 and above)	Calibrator 1	7.5
	Calibrator 2	19.4
	Calibrator 3	44.7
	Calibrator 4	91.0
	Calibrator 5	153.2
	Calibrator 6	177.3