



Table of CRP Control Reference Value

C-Reactive Protein Control

Code: DM3513

 2023112101

 20241120

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DH-100T、DH-102T、DH-103T、 DH-105T、DH-106T、DH-108T、 DH-2000T[H1T1]、 DH-2000T[H2T1]、 DH-2000T[H1T2]、 DH-2000T[H2T2] (Applicable to technical document version A1.0)	Control 1	2.70	1.89-3.51	±0.81
	Control 2	14.32	10.03-18.61	±4.29
	Control 3	40.19	28.13-52.25	±12.06
DH-100T、DH-102T、DH-103T、 DH-105T、DH-106T、DH-108T、 DH-2000T[H1T1]、 DH-2000T[H2T1]、 DH-2000T[H1T2]、 DH-2000T[H2T2] (Applicable to technical document version A2.0)	Control 1	2.94	2.06-3.82	±0.88
	Control 2	14.02	9.81-18.23	±4.21
	Control 3	40.59	28.41-52.77	±12.18