



C-Reactive Protein Control for Vet Reference Value

C-Reactive Protein for Vet Control

Code: DM3702

 2023091201

 20240911

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DF66X Vet、DF68X Vet (Applicable to technical document versions A1.0 and above)	Control 1	16.1	11.3-20.9	±4.8
	Control 2	56.0	39.2-72.8	±16.8
DH10Vet (Applicable to technical document versions A1.0 and above)	Control 1	17.4	12.2-22.6	±5.2
	Control 2	56.1	39.3-72.9	±16.8