

DM-CAL

HEMATOLOGY CALIBRATOR

LOT N0124

2023-12-13

2024-02-10

| Applicable Instruments | Parameter | Unit | Reference Value |
|---|------------|--------------------|-----------------|
| DYMIND DP-H10、 DP-H12、 DP-H13、 DP-H15、 DP-H16、 DP-H17 (Technical File Version A1.10 or higher) ; | WBC | $\times 10^9/L$ | 9.69 |
| | RBC | $\times 10^{12}/L$ | 4.67 |
| | HGB | g/L | 136 |
| | MCV | fL | 90.9 |
| | PLT | $\times 10^9/L$ | 279 |

【NOTE】

1. The calibrator should be stored in refrigerator (2°C~8°C). After opening, it will keep stable for 7 days when it is stored airtight at 2°C~8°C.
2. Before mixing and running the calibrator after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15°C~30°C).
3. Calibrator must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
4. After using, put the calibrator back into the refrigerator to prevent contamination and evaporation.