

	文件类型	技术文件	文件版本	V5.1.0	
	项目代号	DM3113	页 码	第 0 页 共 2 页	
	文件编号	GL DM-56-08-en			
文件名称: <p style="text-align: center;"><b>Table of CRP Control Reference Value</b></p>					
适用范围: <p style="text-align: center;">DM3113 项目</p>					
编制: 曾婷婷 日期: 2024.01.29		审核: 田文飞 日期: 2024.01.29		批准: 田文飞 日期: 2024.01.29	
修 订 记 录					
版本	TCN/ECR/P CN 编号	修订内容概要	修改人	批准人	批准日期
V5.1.0	/	新建	曾婷婷	田文飞	2024.01.29
发至:	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input type="checkbox"/> 营销部 <input type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input type="checkbox"/> 设备研发 <input type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input type="checkbox"/> 试剂质检 <input type="checkbox"/> 设备质 检				
由试剂研发					

## Table of CRP Control Reference Value

### C-Reactive Protein Control

**Code: DM3113**

LOT 2024012201

🕒 20250121

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.9)	<b>Control 1</b>	2.6	1.8-3.4	±0.8
	<b>Control 2</b>	12.7	8.9-16.5	±3.8
	<b>Control 3</b>	28.9	20.2-37.6	±8.7
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above with heating)	<b>Control 1</b>	3.6	2.5-4.7	±1.1
	<b>Control 2</b>	16.9	11.8-22.0	±5.1
	<b>Control 3</b>	37.8	26.5-49.1	±11.3
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above without heating)	<b>Control 1</b>	3.6	2.5-4.7	±1.1
	<b>Control 2</b>	16.8	11.8-21.8	±5.0
	<b>Control 3</b>	37.9	26.5-49.3	±11.4
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A5.11 and below)	<b>Control 1</b>	3.3	2.3-4.3	±1.0
	<b>Control 2</b>	17.5	12.2-22.8	±5.3
	<b>Control 3</b>	36.4	25.5-47.3	±10.9
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above with heating)	<b>Control 1</b>	3.8	2.7-4.9	±1.1
	<b>Control 2</b>	17.4	12.2-22.6	±5.2
	<b>Control 3</b>	40.8	28.6-53.0	±12.2
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above without heating)	<b>Control 1</b>	3.6	2.5-4.7	±1.1
	<b>Control 2</b>	17.8	12.5-23.1	±5.3
	<b>Control 3</b>	40.4	28.3-52.5	±12.1

D6-CRP, DH76CRP (Applicable to technical document version B2.0)	<b>Control 1</b>	3.8	2.7-4.9	$\pm 1.1$
	<b>Control 2</b>	16.6	11.6-21.6	$\pm 5.0$
	<b>Control 3</b>	40.7	28.5-52.9	$\pm 12.2$
D6-CRP, DH76CRP (Applicable to technical document version B2.2 and above)	<b>Control 1</b>	4.1	2.9-5.3	$\pm 1.2$
	<b>Control 2</b>	18.2	12.7-23.7	$\pm 5.5$
	<b>Control 3</b>	39.5	27.6-51.4	$\pm 11.9$
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A8.0 and below with heating)	<b>Control 1</b>	3.9	2.7-5.1	$\pm 1.2$
	<b>Control 2</b>	17.4	12.2-22.6	$\pm 5.2$
	<b>Control 3</b>	41.3	28.9-53.7	$\pm 12.4$
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A8.0 and below without heating)	<b>Control 1</b>	3.7	2.6-4.8	$\pm 1.1$
	<b>Control 2</b>	17.7	12.4-23.0	$\pm 5.3$
	<b>Control 3</b>	42.3	29.6-55.0	$\pm 12.7$
DM79X, DM78X, DM72X (Applicable to technical document version A1.0 and above)	<b>Control 1</b>	4.1	2.9-5.3	$\pm 1.2$
	<b>Control 2</b>	17.9	12.5-23.3	$\pm 5.4$
	<b>Control 3</b>	43.5	30.4-56.6	$\pm 13.1$
DH36X, DH33X, DH31X (Applicable to technical document versions A1.0 and above)	<b>Control 1</b>	4.0	2.8-5.2	$\pm 1.2$
	<b>Control 2</b>	17.5	12.2-22.8	$\pm 5.3$
	<b>Control 3</b>	39.8	27.9-51.7	$\pm 11.9$
DP-H10, DP-H12, DP-H13, DP-H15, DP-H16, DP-H17 (Applicable to technical document versions A1.0 and above)	<b>Control 1</b>	3.9	2.7-5.1	$\pm 1.2$
	<b>Control 2</b>	17.8	12.5-23.1	$\pm 5.3$
	<b>Control 3</b>	36.3	25.4-47.2	$\pm 10.9$