



Table of CRP Calibrator for Vet Reference Value

C-Reactive Protein Calibrator for Vet

Code: DM3702

 2024012401

 20250123

Model	Level (multi-point)	Reference Value (mg/L)
DF66X Vet DF68X Vet (Applicable to technical document versions A1.0 and above)	Calibrator 1	16.3
	Calibrator 2	23.8
	Calibrator 3	56.2
	Calibrator 4	94.1
	Calibrator 5	113.2
	Calibrator 6	155.9
DH10Vet (Applicable to technical document versions A1.0 and above)	Calibrator 1	15.7
	Calibrator 2	26.4
	Calibrator 3	60.8
	Calibrator 4	87.4
	Calibrator 5	119.8
	Calibrator 6	147.8