



## Table of CRP Calibrator Reference Value

### C-Reactive Protein Calibrator

**Code: DM3113**

LOT 2024022001

 20250219

Model	Level (multi-point)	Reference Value (mg/L)
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.9)	Calibrator 1	4.3
	Calibrator 2	20.8
	Calibrator 3	54.4
	Calibrator 4	109.6
	Calibrator 5	140.5
	Calibrator 6	170.3
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above with heating)	Calibrator 1	5.4
	Calibrator 2	20.6
	Calibrator 3	46.9
	Calibrator 4	98.9
	Calibrator 5	139.5
	Calibrator 6	176.4
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above without heating)	Calibrator 1	5.4
	Calibrator 2	21.3
	Calibrator 3	46.7
	Calibrator 4	99.8
	Calibrator 5	135.6
	Calibrator 6	175.9
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A5.11 and below)	Calibrator 1	4.7
	Calibrator 2	21.1
	Calibrator 3	51.3
	Calibrator 4	99.5
	Calibrator 5	134.7
	Calibrator 6	166.2
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above with heating)	Calibrator 1	5.4
	Calibrator 2	20.3
	Calibrator 3	51.2
	Calibrator 4	109.1
	Calibrator 5	143.2
	Calibrator 6	168.1

DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above without heating)	<b>Calibrator 1</b>	5.1
	<b>Calibrator 2</b>	21.6
	<b>Calibrator 3</b>	52.3
	<b>Calibrator 4</b>	102.8
	<b>Calibrator 5</b>	134.7
	<b>Calibrator 6</b>	172.8
D6-CRP, DH76CRP (Applicable to technical document version B2.0)	<b>Calibrator 1</b>	5.1
	<b>Calibrator 2</b>	20.6
	<b>Calibrator 3</b>	54.5
	<b>Calibrator 4</b>	108.9
	<b>Calibrator 5</b>	137.5
	<b>Calibrator 6</b>	178.4
D6-CRP, DH76CRP (Applicable to technical document version B2.2 and above)	<b>Calibrator 1</b>	5.5
	<b>Calibrator 2</b>	21.9
	<b>Calibrator 3</b>	52.3
	<b>Calibrator 4</b>	104.6
	<b>Calibrator 5</b>	140.6
	<b>Calibrator 6</b>	170.5
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A8.1 and below with heating)	<b>Calibrator 1</b>	5.5
	<b>Calibrator 2</b>	22.5
	<b>Calibrator 3</b>	50.1
	<b>Calibrator 4</b>	108.4
	<b>Calibrator 5</b>	144.4
	<b>Calibrator 6</b>	175.3
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A8.1 and below without heating)	<b>Calibrator 1</b>	5.4
	<b>Calibrator 2</b>	22.4
	<b>Calibrator 3</b>	49.6
	<b>Calibrator 4</b>	107.9
	<b>Calibrator 5</b>	140.1
	<b>Calibrator 6</b>	171.4
DM79X, DM78X, DM72X (Applicable to technical document version A1.0 and above)	<b>Calibrator 1</b>	5.9
	<b>Calibrator 2</b>	22.4
	<b>Calibrator 3</b>	54.8
	<b>Calibrator 4</b>	105.2
	<b>Calibrator 5</b>	143.1
	<b>Calibrator 6</b>	170.2

DH36X、DH33X、DH31X (Applicable to technical document versions A1.0 and above)	<b>Calibrator 1</b>	5.3
	<b>Calibrator 2</b>	21.7
	<b>Calibrator 3</b>	52.9
	<b>Calibrator 4</b>	107.6
	<b>Calibrator 5</b>	141.5
	<b>Calibrator 6</b>	176.8
DP-H10、DP-H12、DP-H13、DP-H15、DP-H16、 DP-H17 (Applicable to technical document versions A1.0 and above)	<b>Calibrator 1</b>	6.1
	<b>Calibrator 2</b>	24.4
	<b>Calibrator 3</b>	53.2
	<b>Calibrator 4</b>	97.7
	<b>Calibrator 5</b>	147.0
	<b>Calibrator 6</b>	170.6