

	文件类型	技术文档	文件版本	V3.1.1	
	项目代号	DM 3514	页 码	第 0 页 共 1 页	
	文件编号	GL DM-56-15-en			
文件名称： <p style="text-align: center;"><b>Table of Serum Amyloid A Protein Calibrator Reference Value</b></p>					
适用范围： <p style="text-align: center;">DM 3514 项目</p>					
编制：刘嘉 日期：2024.05.08		审核：刘金盈 日期：2024.05.08		批准：郑剑通 日期：2024.05.08	
<b>修 订 记 录</b>					
版本	TCN/ECR/PCN 编号	修订内容概要	修改人	批准人	批准日期
V3.1.0	/	新增	秦俊飞	田文飞	2024.04.15
V3.1.1	/	增加1117仪器参考值	刘嘉	郑剑通	2024.05.08
发至：	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input checked="" type="checkbox"/> 营销部 <input checked="" type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input checked="" type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input checked="" type="checkbox"/> 设备研发 <input checked="" type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input checked="" type="checkbox"/> 设备生产 <input checked="" type="checkbox"/> 试剂生产 <input checked="" type="checkbox"/> 试剂质检 <input checked="" type="checkbox"/> 设备质检				
由设转					

## Table of Serum Amyloid A Protein Calibrator Reference Value

### Serum Amyloid A Protein Calibrator

**Code: DM3514**

**LOT** 2024032701

 20250326

Model	Level	Reference Value (mg/L)
DH-100T、DH-102T、DH-103T、DH-105T、DH-106T、 DH-108T、DH-2000T[H1T1]、DH-2000T[H2T1]、 DH-2000T[H1T2]、DH-2000T[H2T2] (Applicable to technical document version A1.0 and above)	<b>Calibrator 1</b>	9.89
	<b>Calibrator 2</b>	21.40
	<b>Calibrator 3</b>	51.66
	<b>Calibrator 4</b>	105.44
	<b>Calibrator 5</b>	185.55
	<b>Calibrator 6</b>	236.98
	<b>Calibrator 7</b>	275.49
DH-800 CRP series DH-800 CS series (Technical File Version A1.3 or higher)	<b>Calibrator 1</b>	9.20
	<b>Calibrator 2</b>	18.00
	<b>Calibrator 3</b>	49.20
	<b>Calibrator 4</b>	101.00
	<b>Calibrator 5</b>	179.00
	<b>Calibrator 6</b>	243.50
	<b>Calibrator 7</b>	251.00