


Table of Serum Amyloid A Protein Control Reference Value

Serum Amyloid A Protein Control

Code: DM3514

LOT 2024051801

 20250517

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DH-100T、DH-102T、DH-103T、 DH-105T、DH-106T、DH-108T、 DH-2000T[H1T1]、 DH-2000T[H2T1]、 DH-2000T[H1T2]、 DH-2000T[H2T2] (Applicable to technical document version A1.0 and above)	Control 1	7.34	4.64-10.04	±2.70
	Control 2	51.74	36.80-66.68	±14.94
DH-800 CRP series DH-800 CS series (Technical File Version A1.3 or higher)	Control 1	6.58	4.61~8.56	±1.97
	Control 2	49.60	34.72~64.48	±14.88