

	文件类型	技术文件	文件版本	V5.5.0	
	项目代号	DM3113	页 码	第 0 页 共 2 页	
	文件编号	GL DM-56-08-en			
文件名称: <p style="text-align: center;">Table of CRP Control Reference Value</p>					
适用范围: <p style="text-align: center;">DM3113 项目</p>					
编制: 黄子苾 日期: 2024.07.26		审核: 曾婷婷 日期: 2024.07.26		批准: 田文飞 日期: 2024.07.26	
修 订 记 录					
版本	TCN/ECR/P CN 编号	修订内容概要	修改人	批准人	批准日期
V5.5.0	/	新建	黄子苾	田文飞	2024.07.26
发至:	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input type="checkbox"/> 营销部 <input type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input type="checkbox"/> 设备研发 <input type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input type="checkbox"/> 试剂质检 <input type="checkbox"/> 设备质 检				
由试剂研发					

Table of CRP Control Reference Value

C-Reactive Protein Control

Code: DM3113

LOT 2024070201

🕒 20250701

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.9)	Control 1	2.4	1.7-3.1	±0.7
	Control 2	14.1	9.9-18.3	±4.2
	Control 3	35.9	25.1-46.7	±10.8
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above with heating)	Control 1	2.9	2.0-3.8	±0.9
	Control 2	13.4	9.4-17.4	±4.0
	Control 3	34.5	24.2-44.9	±10.4
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above without heating)	Control 1	2.8	2.0-3.6	±0.8
	Control 2	13.7	9.6-17.8	±4.1
	Control 3	36.9	25.8-48.0	±11.1
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A5.11 and below)	Control 1	1.9	1.3-2.5	±0.6
	Control 2	12.8	9.0-16.6	±3.8
	Control 3	36.2	25.3-47.1	±10.9
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above with heating)	Control 1	3.0	2.1-3.9	±0.9
	Control 2	15.2	10.6-19.8	±4.6
	Control 3	38.9	27.2-50.6	±11.7
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above without heating)	Control 1	3.5	2.5-4.6	±1.1
	Control 2	14.8	10.4-19.2	±4.4
	Control 3	38.5	27.0-50.1	±11.6
D6-CRP, DH76CRP	Control 1	2.7	1.9-3.5	±0.8

(Applicable to technical document version B2.0)	Control 2	13.4	9.4-17.4	±4.0
	Control 3	37.4	26.2-48.6	±11.2
D6-CRP, DH76CRP (Applicable to technical document version B2.2 and above)	Control 1	2.3	1.6-3.0	±0.7
	Control 2	14.5	10.2-18.9	±4.4
	Control 3	35.6	24.92-46.28	±10.7
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A8.1 and below with heating)	Control 1	3.6	2.52-4.68	±1.1
	Control 2	14.3	10.0-18.6	±4.3
	Control 3	37.9	26.5-49.3	±11.4
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A8.1 and below without heating)	Control 1	2.8	2.0-3.6	±0.8
	Control 2	13.6	9.5-17.7	±4.1
	Control 3	47.1	33.0-61.2	±14.1
DM79X, DM78X, DM72X (Applicable to technical document version A1.0 and above)	Control 1	3.7	2.6-4.8	±1.1
	Control 2	16.7	11.7-21.7	±5.0
	Control 3	40.4	28.3-52.5	±12.1
DH36X, DH33X, DH31X (Applicable to technical document versions A1.0 and above)	Control 1	2.8	2.0-3.6	±0.8
	Control 2	13.7	9.6-17.8	±4.1
	Control 3	36.8	25.8-47.8	±11.0
DP-H10, DP-H12, DP-H13, DP-H15, DP-H16, DP-H17 (Applicable to technical document versions A1.0 and above)	Control 1	4.1	2.87-5.33	±1.2
	Control 2	17.2	12.0-22.4	±5.2
	Control 3	42.2	29.5-54.9	±12.7
DH-800 CRP series DH-800 CS series (Technical File Version A1.3 or higher)	Control 1	3.1	2.2-4.0	±0.9
	Control 2	15.0	10.5-19.5	±4.5
	Control 3	33.4	23.4-43.4	±10.0