
	文件类型	技术文件	文件版本	V5.5.0	
	项目代号	DM3113	页 码	第 0 页 共 3 页	
	文件编号	GL DM-56-09-en			
文件名称: <p style="text-align: center;">Table of CRP Calibrator Reference Value</p>					
适用范围: <p style="text-align: center;">DM3113 项目</p>					
编制: 黄子苾 日期: 2024.07.26		审核: 曾婷婷 日期: 2024.07.26		批准: 魏保立 日期: 2024.07.26	
修 订 记 录					
版本	TCN/ECR/P CN 编号	修订内容概要	修改人	批准人	批准日期
V5.5.0	/	新建	黄子苾	田文飞	2024.07.26
发至:	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input type="checkbox"/> 营销部 <input type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input type="checkbox"/> 设备研发 <input type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input type="checkbox"/> 试剂质检 <input type="checkbox"/> 设备质 检				
由试剂研发					

Table of CRP Calibrator Reference Value

C-Reactive Protein Calibrator

Code: DM3113

 2024070201

 20250701

Model	Level (multi-point)	Reference Value (mg/L)
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.9)	Calibrator 1	4.6
	Calibrator 2	18.5
	Calibrator 3	42.1
	Calibrator 4	86.6
	Calibrator 5	110.7
	Calibrator 6	138.6
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above with heating)	Calibrator 1	5.0
	Calibrator 2	17.4
	Calibrator 3	46.1
	Calibrator 4	96.6
	Calibrator 5	133.6
	Calibrator 6	184.3
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above without heating)	Calibrator 1	5.3
	Calibrator 2	19.5
	Calibrator 3	47.2
	Calibrator 4	100.6
	Calibrator 5	152.3
	Calibrator 6	194.5
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A5.11 and below)	Calibrator 1	4.2
	Calibrator 2	17.8
	Calibrator 3	46.8
	Calibrator 4	101.3
	Calibrator 5	133.2
	Calibrator 6	196.8
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above with heating)	Calibrator 1	5.2
	Calibrator 2	21.2
	Calibrator 3	50.6
	Calibrator 4	97.5
	Calibrator 5	147.8
	Calibrator 6	181.2

DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above without heating)	Calibrator 1	5.6
	Calibrator 2	18.8
	Calibrator 3	51.8
	Calibrator 4	99.6
	Calibrator 5	148.3
	Calibrator 6	190.3
D6-CRP, DH76CRP (Applicable to technical document version B2.0)	Calibrator 1	5.2
	Calibrator 2	19.5
	Calibrator 3	48.1
	Calibrator 4	105.1
	Calibrator 5	151.3
	Calibrator 6	182.3
D6-CRP, DH76CRP (Applicable to technical document version B2.2 and above)	Calibrator 1	5.4
	Calibrator 2	18.9
	Calibrator 3	47.1
	Calibrator 4	97.2
	Calibrator 5	140.2
	Calibrator 6	174.6
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A8.1 and below with heating)	Calibrator 1	5.3
	Calibrator 2	19.3
	Calibrator 3	53.6
	Calibrator 4	105.9
	Calibrator 5	156.7
	Calibrator 6	192.3
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A8.1 and below without heating)	Calibrator 1	5.5
	Calibrator 2	22.8
	Calibrator 3	51.8
	Calibrator 4	108.6
	Calibrator 5	148.7
	Calibrator 6	195.6
DM79X, DM78X, DM72X (Applicable to technical document version A1.0 and above)	Calibrator 1	6.4
	Calibrator 2	21.6
	Calibrator 3	51.8
	Calibrator 4	103.9
	Calibrator 5	141.2
	Calibrator 6	182.6

DH36X、DH33X、DH31X (Applicable to technical document versions A1.0 and above)	Calibrator 1	5.1
	Calibrator 2	19.4
	Calibrator 3	50.6
	Calibrator 4	105.7
	Calibrator 5	145.3
	Calibrator 6	179.3
DP-H10、DP-H12、DP-H13、DP-H15、DP-H16、 DP-H17 (Applicable to technical document versions A1.0 and above)	Calibrator 1	6.4
	Calibrator 2	21.8
	Calibrator 3	51.2
	Calibrator 4	102.3
	Calibrator 5	163.1
	Calibrator 6	196.6
DH-800 CRP 系列 DH-800 CS 系列 (Applicable to technical document versions A1.3 and above)	Calibrator 1	5.5
	Calibrator 2	20.1
	Calibrator 3	41.1
	Calibrator 4	85.2
	Calibrator 5	112.0
	Calibrator 6	154.2