




## Table of CRP Calibrator Reference Value

### C-Reactive Protein Calibrator

**Code: DM3113**

 2024070201

 20250701

Model	Level (multi-point)	Reference Value (mg/L)
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.9)	Calibrator 1	4.6
	Calibrator 2	18.5
	Calibrator 3	42.1
	Calibrator 4	86.6
	Calibrator 5	110.7
	Calibrator 6	138.6
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above with heating)	Calibrator 1	5.0
	Calibrator 2	17.4
	Calibrator 3	46.1
	Calibrator 4	96.6
	Calibrator 5	133.6
	Calibrator 6	184.3
DH51CRP, DH53CRP, DH56CRP D1-CRP, D3-CRP, D5-CRP (Applicable to technical document version A11.10 and above without heating)	Calibrator 1	5.3
	Calibrator 2	19.5
	Calibrator 3	47.2
	Calibrator 4	100.6
	Calibrator 5	152.3
	Calibrator 6	194.5
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A5.11 and below)	Calibrator 1	4.2
	Calibrator 2	17.8
	Calibrator 3	46.8
	Calibrator 4	101.3
	Calibrator 5	133.2
	Calibrator 6	196.8
DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above with heating)	Calibrator 1	5.2
	Calibrator 2	21.2
	Calibrator 3	50.6
	Calibrator 4	97.5
	Calibrator 5	147.8
	Calibrator 6	181.2

DH71CRP, DH73CRP D2-CRP, D7-CRP (Applicable to technical document versions A6.0 and above without heating)	<b>Calibrator 1</b>	5.6
	<b>Calibrator 2</b>	18.8
	<b>Calibrator 3</b>	51.8
	<b>Calibrator 4</b>	99.6
	<b>Calibrator 5</b>	148.3
	<b>Calibrator 6</b>	190.3
D6-CRP, DH76CRP (Applicable to technical document version B2.0)	<b>Calibrator 1</b>	5.2
	<b>Calibrator 2</b>	19.5
	<b>Calibrator 3</b>	48.1
	<b>Calibrator 4</b>	105.1
	<b>Calibrator 5</b>	151.3
	<b>Calibrator 6</b>	182.3
D6-CRP, DH76CRP (Applicable to technical document version B2.2 and above)	<b>Calibrator 1</b>	5.4
	<b>Calibrator 2</b>	18.9
	<b>Calibrator 3</b>	47.1
	<b>Calibrator 4</b>	97.2
	<b>Calibrator 5</b>	140.2
	<b>Calibrator 6</b>	174.6
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A8.1 and below with heating)	<b>Calibrator 1</b>	5.3
	<b>Calibrator 2</b>	19.3
	<b>Calibrator 3</b>	53.6
	<b>Calibrator 4</b>	105.9
	<b>Calibrator 5</b>	156.7
	<b>Calibrator 6</b>	192.3
DF1-CRP, DF3-CRP, DF5-CRP DF50CRP, DF52CRP, DF53CRP (Applicable to technical document versions A8.1 and below without heating)	<b>Calibrator 1</b>	5.5
	<b>Calibrator 2</b>	22.8
	<b>Calibrator 3</b>	51.8
	<b>Calibrator 4</b>	108.6
	<b>Calibrator 5</b>	148.7
	<b>Calibrator 6</b>	195.6
DM79X, DM78X, DM72X (Applicable to technical document version A1.0 and above)	<b>Calibrator 1</b>	6.4
	<b>Calibrator 2</b>	21.6
	<b>Calibrator 3</b>	51.8
	<b>Calibrator 4</b>	103.9
	<b>Calibrator 5</b>	141.2
	<b>Calibrator 6</b>	182.6

DH36X、DH33X、DH31X (Applicable to technical document versions A1.0 and above)	<b>Calibrator 1</b>	5.1
	<b>Calibrator 2</b>	19.4
	<b>Calibrator 3</b>	50.6
	<b>Calibrator 4</b>	105.7
	<b>Calibrator 5</b>	145.3
	<b>Calibrator 6</b>	179.3
DP-H10、DP-H12、DP-H13、DP-H15、DP-H16、 DP-H17 (Applicable to technical document versions A1.0 and above)	<b>Calibrator 1</b>	6.4
	<b>Calibrator 2</b>	21.8
	<b>Calibrator 3</b>	51.2
	<b>Calibrator 4</b>	102.3
	<b>Calibrator 5</b>	163.1
	<b>Calibrator 6</b>	196.6
DH-800 CRP series DH-800 CS series (Technical File Version A1.3 or higher)	<b>Calibrator 1</b>	5.5
	<b>Calibrator 2</b>	20.1
	<b>Calibrator 3</b>	41.1
	<b>Calibrator 4</b>	85.2
	<b>Calibrator 5</b>	112.0
	<b>Calibrator 6</b>	154.2