



Table of CRP Control Reference Value

C-Reactive Protein Control

Code: DM3513

 2024032101

 20250320

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DH-100T、DH-102T、DH-103T、 DH-105T、DH-106T、DH-108T、 DH-2000T[H1T1]、 DH-2000T[H2T1]、 DH-2000T[H1T2]、 DH-2000T[H2T2] (Applicable to technical document version A1.0)	Control 1	3.27	2.29-4.25	±0.98
	Control 2	17.90	12.53-23.27	±5.37
	Control 3	44.63	31.24-58.02	±13.39
DH-100T、DH-102T、DH-103T、 DH-105T、DH-106T、DH-108T、 DH-2000T[H1T1]、 DH-2000T[H2T1]、 DH-2000T[H1T2]、 DH-2000T[H2T2] (Applicable to technical document version A2.0 and above)	Control 1	3.56	2.49-4.63	±1.07
	Control 2	17.39	12.17-22.61	±5.22
	Control 3	44.88	31.42-58.34	±13.46