


Table of CRP Control Reference Value

C-Reactive Protein Control

Code: DM3513

LOT 2024070201

 20250701

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DH-100T、DH-102T、DH-103T、 DH-105T、DH-106T、DH-108T、 DH-2000T[H1T1]、 DH-2000T[H2T1]、 DH-2000T[H1T2]、 DH-2000T[H2T2] (Applicable to technical document version A1.0)	Control 1	3.80	2.66-4.94	±1.14
	Control 2	19.60	13.72-25.48	±5.88
	Control 3	41.50	29.05-53.95	±12.45
DH-100T、DH-102T、DH-103T、 DH-105T、DH-106T、DH-108T、 DH-2000T[H1T1]、 DH-2000T[H2T1]、 DH-2000T[H1T2]、 DH-2000T[H2T2] (Applicable to technical document version A2.0 and above)	Control 1	3.80	2.66-4.94	±1.14
	Control 2	16.90	11.83-21.97	±5.07
	Control 3	39.60	27.72-51.48	±11.88