





## Table of Serum Amyloid A Protein Control Reference Value

### Serum Amyloid A Protein Control

**Code: DM3514**

 2024032701

 20250326

Model	Level	Reference Value (mg/L)	Reference Range (mg/L)	Deviation Range (mg/L)
DH-100T、DH-102T、DH-103T、 DH-105T、DH-106T、DH-108T、 DH-2000T[H1T1]、 DH-2000T[H2T1]、 DH-2000T[H1T2]、 DH-2000T[H2T2] (Applicable to technical document version A1.0 and above)	<b>Control 1</b>	8.00	5.60-10.40	±2.40
	<b>Control 2</b>	52.46	36.72-68.20	±15.74
DH-800 CRP series DH-800 CS series ( Technical File Version A1.3 or higher )	<b>Control 1</b>	7.80	5.46-10.14	±2.34
	<b>Control 2</b>	48.76	34.13-63.38	±14.63