





# Fibrinogen Calibrator Reference Value Table

## Fibrinogen Calibrator

### CALIBRATOR

LOT 2024120501

 20241205 20261204

Applicable Instrument	Unit	Principle	Reference Value
CA1200 CA1203 CA1205 DC1210	g/L	Optical method	3.11
DC1250 DC1280	g/L	Magnetic bead	3.06

**【NOTE】**

1. The calibrator should be stored in the refrigerator at 2 °C - 8 °C. After reconstitution, the calibrator can be stable for 4 hours at room temperature (10 °C - 30 °C), and can be stable for 8 hours at 2 °C - 8 °C.
2. The calibrator cannot be frozen after reconstitution
3. The calibrator is for IVD only.