


Table of C-Reactive Protein Calibrator (CRP) Reference Value

C-Reactive Protein Calibrator

LOT 2024111501

 2025-11-14

Applicable Instruments	DH51CRP, DH53CRP, DH56CRP, D1-CRP, D3-CRP, D5-CRP (Technical file version A11.10 or higher with heating)	DH51CRP, DH53CRP, DH56CRP, D1-CRP, D3-CRP, D5-CRP (Technical file version A11.10 or higher without heating)	DH71CRP, DH73CRP, D2-CRP, D7-CRP (Technical file version A6.0 or higher with heating)	DH71CRP, DH73CRP, D2-CRP, D7-CRP (Technical file version A6.0 or higher without heating)	D6-CRP, DH76CRP (Technical file version B2.0)	D6-CRP, DH76CRP (Technical file version B2.2 or higher)	DF1-CRP, DF3-CRP, DF5-CRP, DF50CRP, DF52CRP, DF53CRP (Technical file version A6.0 or higher with heating)
Level (multi-point)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)
Calibrator 1	5.3	5.3	4.9	5.2	4.9	4.7	5.0
Calibrator 2	17.9	18.9	19.6	18.8	18.6	18.0	17.6
Calibrator 3	50.5	55.9	55.0	54.1	50.7	54.1	51.5
Calibrator 4	97.9	98.0	94.6	97.9	92.3	86.1	96.7
Calibrator 5	159.8	168.4	162.0	158.7	145.1	131.6	167.0
Calibrator 6	189.7	196.6	198.0	190.0	173.5	178.4	191.8
Applicable Instruments	DF1-CRP, DF3-CRP, DF5-CRP, DF50CRP, DF52CRP, DF53CRP (Technical file versions A8.1 and below without heating)	DM79X, DM78X, DM72X (Technical file version A1.0 or higher)	DH36X, DH33X, DH31X (Technical file versions A1.0 or higher)	DP-H10, DP-H12, DP-H13, DP-H15, DP-H16, DP-H17 (Technical file versions A1.0 or higher)	DH-800 CRP series DH-800 CS series (Technical file Version A1.3 or higher)	DH-2000T series , DH-100T, DH-102T DH-103T, DH-105T DH-106T, DH-108T (Technical file version A2.0 or higher)	/
Level (multi-point)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	/
Calibrator 1	5.8	5.0	5.1	6.6	5.01	4.71	/
Calibrator 2	20.3	18.7	18.2	25.0	16.82	18.54	/
Calibrator 3	54.5	46.9	55.2	50.8	44.20	49.32	/
Calibrator 4	108.0	87.6	96.0	111.4	77.24	93.37	/
Calibrator 5	161.4	151.9	167.1	156.2	115.45	142.66	/
Calibrator 6	177.4	184.4	192.0	191.1	145.73	178.59	/