


Table of C-Reactive Protein Calibrator (CRP) Reference Value

C-Reactive Protein Calibrator

LOT 2024111501

 2025-11-14

Applicable Instruments	DH51CRP, DH53CRP DH56CRP, D1-CRP D3-CRP, D5-CRP (Technical file version A11.9)	/	/	/	/	/	/
Level (multi-point)	Reference Value (mg/L)	/	/	/	/	/	/
Calibrator 1	4.4	/	/	/	/	/	/
Calibrator 2	10.7	/	/	/	/	/	/
Calibrator 3	15.5	/	/	/	/	/	/
Calibrator 4	41.0	/	/	/	/	/	/
Calibrator 5	85.0	/	/	/	/	/	/
Calibrator 6	145.9	/	/	/	/	/	/