
	文件类型	技术文件	文件版本	V5.9.0		
	项目代号	/	页码	第0页 共1页		
	文件编号	GL DM-56-07-en				
文件名称： <p style="text-align: center;">Table of C-Reactive Protein Calibrator (CRP) Reference Value 2024123001</p> 适用范围： DM3113项目、DM3513项目						
编制：刘嘉 日期：2025.03.07		审核：刘金盈 日期：2025.03.07		批准：程军 日期：2025.03.07		
修 订 记 录						
版本	TCN/ECR/PCN编号	修订内容概要		修改人	批准人	批准日期
V5.9.0	/	新建		刘嘉	程军	2025.03.07
发至：	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input checked="" type="checkbox"/> 营销部 <input checked="" type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input checked="" type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input checked="" type="checkbox"/> 设备研发 <input checked="" type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input checked="" type="checkbox"/> 试剂质检 <input checked="" type="checkbox"/> 设备质检					
由试剂研发						

Table of C-Reactive Protein Calibrator (CRP) Reference Value

C-Reactive Protein Calibrator

LOT 2024123001

 2025-12-29

Applicable Instruments	DH51CRP, DH53CRP, DH56CRP, D1-CRP, D3-CRP, D5-CRP (Technical file version A11.10 or higher with heating)	DH51CRP, DH53CRP, DH56CRP, D1-CRP, D3-CRP, D5-CRP (Technical file version A11.10 or higher without heating)	DH71CRP, DH73CRP, D2-CRP, D7-CRP (Technical file version A6.0 or higher with heating)	DH71CRP, DH73CRP, D2-CRP, D7-CRP (Technical file version A6.0 or higher without heating)	D6-CRP, DH76CRP (Technical file version B2.0)	D6-CRP, DH76CRP (Technical file version B2.2 or higher)	DF1-CRP, DF3-CRP, DF5-CRP, DF50CRP, DF52CRP, DF53CRP (Technical file version A6.0 or higher with heating)
Level (multi-point)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)
Calibrator 1	7.9	7.8	5.0	6.6	5.9	6.3	7.7
Calibrator 2	22.9	24.9	20.2	22.9	21.5	22.2	24.9
Calibrator 3	60.8	60.8	53.5	58.5	53.6	53.1	65.1
Calibrator 4	140.9	135.2	115.1	127.5	121.5	111.0	164.5
Calibrator 5	158.5	168.7	159.1	171.7	169.3	171.4	195.4
Calibrator 6	188.0	208.7	216.7	208.1	204.3	198.9	212.5
Applicable Instruments	DF1-CRP, DF3-CRP, DF5-CRP, DF50CRP, DF52CRP, DF53CRP (Technical file versions A8.1 and below without heating)	DM79X, DM78X, DM72X (Technical file version A1.0 or higher)	DH36X, DH33X, DH31X (Technical file versions A1.0 or higher)	DP-H10, DP-H12, DP-H13, DP-H15, DP-H16, DP-H17 (Technical file versions A1.0 or higher)	DH-800 CRP series DH-800 CS series (Technical file Version A1.3 or higher)	DH-2000T series , DH-100T, DH-102T DH-103T, DH-105T DH-106T, DH-108T (Technical file version A2.0 or higher)	/
Level (multi-point)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	/
Calibrator 1	7.7	6.00	6.2	6.2	5.65	5.51	/
Calibrator 2	24.5	21.50	21.9	19.1	18.20	20.95	/
Calibrator 3	65.2	50.89	55.1	49.4	47.50	52.55	/
Calibrator 4	139.9	114.90	134.6	95.2	90.80	105.21	/
Calibrator 5	163.1	169.17	170.7	161.7	135.96	159.51	/
Calibrator 6	193.2	202.50	202.2	180.2	166.50	208.57	/