


Fibrin and Fibrinogen Degradation Product Calibrator


Reference Value Table

FDP Calibrator

CALIBRATOR

LOT 2024102901

 20241029

 20261028

Applicable Instrument	Unit	Reference Value
CA1200	μg/mL	FDP Calibrator 1: 4.75
CA1203		FDP Calibrator 2: 20.64
CA1205		FDP Calibrator 3: 40.39
DC1210		FDP Calibrator 4: 73.82
DC1250		FDP Calibrator 5: 107.75
DC1280		

【NOTE】

1. The calibrator should be stored in the refrigerator at 2°C ~ 8°C. After reconstitution, the calibrator can be stable for 8 hours at room temperature (10°C~ 30°C), and can be stable for 12 hours at 2°C ~ 8°C.
2. The calibrator cannot be frozen after reconstitution.
3. The calibrator is for IVD only.