



	文件类型	技术文件	文件版本	V5.3.0		
	项目代号	/	页码	第0页 共2页		
	文件编号	GL DM-56-17-en				
文件名称： <p style="text-align: center;">CONTROL DM-RET RET0525 2025-07-10</p> 适用范围： 1116项目、7003项目						
编制：郭罕星 日期：2025.04.30		审核：刘金盈 日期：2025.04.30		批准：程军 日期：2025.04.30		
修 订 记 录						
版本	TCN/ECR/PCN编号	修订内容概要		修改人	批准人	批准日期
V5.3.0	/	新建		郭罕星	程军	2025.04.30
发至：	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input checked="" type="checkbox"/> 营销部 <input checked="" type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input checked="" type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input checked="" type="checkbox"/> 设备研发 <input checked="" type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input checked="" type="checkbox"/> 试剂质检 <input checked="" type="checkbox"/> 设备质检					
由工程						

DM-RET

HEMATOLOGY CONTROL

Reference Values provided by DYMIND
CONTROL
 2025-04-12

 2025-07-10

Applicable Instruments	Parameter	Unit	LOT RET0525L	LOT RET0525N	LOT RET0525H
DYMIND DH-600 DH-602 DH-605 DH-610 DH-612 DH-615 (Technical File Version A2.0 or higher)	RET#	$\times 10^{12}/L$	0.0416 ± 0.0208	0.2303 ± 0.1316	0.2276 ± 0.1138
	RET%	%	0.76 ± 0.38	4.94 ± 1.98	10.46 ± 4.18
	LFR	%	89.9 ± 27.0	89.9 ± 27.0	89.4 ± 27.0
	MFR	%	8.1 ± 12.0	8.1 ± 12.0	8.6 ± 13.0
	HFR	%	2.0 ± 3.0	2.0 ± 3.0	2.0 ± 3.0
	IRF	%	10.1 ± 15.0	10.1 ± 15.0	10.6 ± 16.0
	RHE	pg	16.3 ± 7.0	16.1 ± 6.0	15.9 ± 6.0

【NOTE】

1. The controls should be stored in refrigerator (2°C~8°C). After opening, it will keep stable for 14 days when it is stored airtight at 2°C~8°C.
2. Before mixing and running the control after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15°C~30°C).
3. Controls must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
4. After using, put the controls back into the refrigerator to prevent contamination and evaporation.



RET0525L




RET0525N




RET0525H

DM-RET

HEMATOLOGY CONTROL

Reference Values provided by DYMIND
CONTROL
 2025-04-12

 2025-07-10

Applicable Instruments	Parameter	Unit	LOT RET0525L	LOT RET0525N	LOT RET0525H
DYMIND DM60 Vet DM61 Vet DM62 Vet DM63 Vet (Technical File Version A1.4 or higher)	RET#	$\times 10^{12}/L$	0.0291 ±0.0233	0.1458 ±0.0729	0.1577 ±0.0789
	RET%	%	0.52 ±0.42	3.21 ±1.61	7.10 ±3.55

【NOTE】

1. The controls should be stored in refrigerator (2°C~8°C). After opening, it will keep stable for 14 days when it is stored airtight at 2°C~8°C.
2. Before mixing and running the control after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15°C~30°C).
3. Controls must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
4. After using, put the controls back into the refrigerator to prevent contamination and evaporation.



RET0525L



RET0525N



RET0525H