




# Table of C-Reactive Protein Calibrator (CRP) Reference Value

## C-Reactive Protein Calibrator

LOT 2025041701

 2026-04-16

Applicable Instruments	DH51CRP, DH53CRP, DH56CRP, D1-CRP, D3-CRP, D5-CRP (Technical file version A11.10 or higher with heating)	DH51CRP, DH53CRP, DH56CRP, D1-CRP, D3-CRP, D5-CRP (Technical file version A11.10 or higher without heating)	DH71CRP, DH73CRP, D2-CRP, D7-CRP (Technical file version A6.0 or higher with heating)	DH71CRP, DH73CRP, D2-CRP, D7-CRP (Technical file version A6.0 or higher without heating)	D6-CRP, DH76CRP (Technical file version B2.0)	D6-CRP, DH76CRP (Technical file version B2.2 or higher)	DF1-CRP, DF3-CRP, DF5-CRP, DF50CRP, DF52CRP, DF53CRP (Technical file version A6.0 or higher with heating)
Level (multi-point)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)
Calibrator 1	4.8	5.2	5.2	4.8	5.2	4.6	4.6
Calibrator 2	16.5	18.3	18.4	17.8	19.0	15.9	15.7
Calibrator 3	42.5	48.2	46.9	43.3	50.3	42.1	43.2
Calibrator 4	94.0	109.8	103.5	92.6	116.6	79.8	100.1
Calibrator 5	126.0	127.0	133.0	113.5	136.5	99.6	148.0
Calibrator 6	146.1	148.0	163.0	157.0	174.1	133.3	164.0
Applicable Instruments	DF1-CRP, DF3-CRP, DF5-CRP, DF50CRP, DF52CRP, DF53CRP (Technical file versions A8.1 and below without heating)	DM79X, DM78X, DM72X (Technical file version A1.0 or higher)	DH36X, DH33X, DH31X (Technical file versions A1.0 or higher)	DP-H10, DP-H12, DP-H13, DP-H15, DP-H16, DP-H17 (Technical file versions A1.0 or higher)	DH-800 CRP series, DH-800 CS series (Technical file Version A1.3 or higher)	DH-2000T series, DH-100T, DH-102T, DH-103T, DH-105T, DH-106T, DH-108T (Technical file version A2.0 or higher)	/
Level (multi-point)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	/
Calibrator 1	5.2	5.4	4.7	6.5	5.38	5.20	
Calibrator 2	17.3	18.6	17.1	21.2	18.94	18.30	
Calibrator 3	46.4	46.4	43.0	46.1	47.34	45.40	
Calibrator 4	95.3	97.5	92.5	101.4	94.60	91.00	
Calibrator 5	138.0	122.1	121.4	126.1	118.60	119.00	
Calibrator 6	158.0	163.6	150.7	144.3	154.00	162.00	