
	文件类型	技术文件	文件版本	V6.1.1		
	项目代号	/	页码	第0页 共1页		
	文件编号	GL DM-56-07-en				
文件名称： <p style="text-align: center;">Table of C-Reactive Protein Calibrator (CRP) Reference Value</p> 适用范围： DM3113项目、DM3107项目、DM3513项目						
编制：郭罕星 日期：2025.05.15		审核：刘金盈 日期：2025.05.15		批准：程军 日期：2025.05.15		
修 订 记 录						
版本	TCN/ECR/PCN编号	修订内容概要		修改人	批准人	批准日期
V6.1.0	/	新建		郭罕星	程军	2025.04.30
V6.1.1	/	更改P101参考值		郭罕星	程军	2025.05.15
发至：	<input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input checked="" type="checkbox"/> 营销部 <input checked="" type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input checked="" type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input checked="" type="checkbox"/> 设备研发 <input checked="" type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input checked="" type="checkbox"/> 试剂质检 <input checked="" type="checkbox"/> 设备质检					
由工程						

Table of C-Reactive Protein Calibrator (CRP) Reference Value

C-Reactive Protein Calibrator

LOT 2025041701

 2026-04-16

Applicable Instruments	DH51CRP, DH53CRP, DH56CRP, D1-CRP, D3-CRP, D5-CRP (Technical file version A11.10 or higher with heating)	DH51CRP, DH53CRP, DH56CRP, D1-CRP, D3-CRP, D5-CRP (Technical file version A11.10 or higher without heating)	DH71CRP, DH73CRP, D2-CRP, D7-CRP (Technical file version A6.0 or higher with heating)	DH71CRP, DH73CRP, D2-CRP, D7-CRP (Technical file version A6.0 or higher without heating)	D6-CRP, DH76CRP (Technical file version B2.0)	D6-CRP, DH76CRP (Technical file version B2.2 or higher)	DF1-CRP, DF3-CRP, DF5-CRP, DF50CRP, DF52CRP, DF53CRP (Technical file version A6.0 or higher with heating)
Level (multi-point)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)
Calibrator 1	4.8	5.2	5.2	4.8	5.2	4.6	4.6
Calibrator 2	16.5	18.3	18.4	17.8	19.0	15.9	15.7
Calibrator 3	42.5	48.2	46.9	43.3	50.3	42.1	43.2
Calibrator 4	94.0	109.8	103.5	92.6	116.6	79.8	100.1
Calibrator 5	126.0	127.0	133.0	113.5	136.5	99.6	148.0
Calibrator 6	146.1	148.0	163.0	157.0	174.1	133.3	164.0
Applicable Instruments	DF1-CRP, DF3-CRP, DF5-CRP, DF50CRP, DF52CRP, DF53CRP (Technical file versions A8.1 and below without heating)	DM79X, DM78X, DM72X (Technical file version A1.0 or higher)	DH36X, DH33X, DH31X (Technical file versions A1.0 or higher)	DP-H10, DP-H12, DP-H13, DP-H15, DP-H16, DP-H17 (Technical file versions A1.0 or higher)	DH-800 CRP series DH-800 CS series (Technical file Version A1.3 or higher)	DH-2000T series , DH-100T, DH-102T DH-103T, DH-105T DH-106T, DH-108T (Technical file version A2.0 or higher)	/
Level (multi-point)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	/
Calibrator 1	5.2	5.40	4.7	6.6	5.38	5.20	
Calibrator 2	17.3	18.60	17.1	21.6	18.94	18.30	
Calibrator 3	46.4	46.36	43.0	46.4	47.34	45.40	
Calibrator 4	95.3	97.50	92.5	101.1	94.60	91.00	
Calibrator 5	138.0	122.10	121.4	141.9	118.60	119.00	
Calibrator 6	158.0	163.59	150.7	181.9	154.00	162.00	