


Table of Serum Amyloid A Protein Control (SAA) Reference Value

Serum Amyloid A Protein Control

LOT 2025050601

 2026-05-05

Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)	Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)
DH-2000T series DH-100T、DH-102T DH-103T、DH-105T DH-106T、DH-108T (Technical File version A2.0 or higher)	Control 1	8.38	±2.51	DH-800 CS series (Technical File Version A1.3 or higher)	Control 1	8.10	±2.43
	Control 2	49.05	±14.72		Control 2	58.25	±17.48