



D-dimer Control Reference Value Table

D-dimer Control

CONTROL

LOT 2025052101

 20250521

 20270520

Applicable Instruments	Unit	Concentration	Reference Value	Deviation Limit
CA1200 CA1203 CA1205	μg/mL (FEU)	Low value	0.40	± 0.12
DC1210 DC1250 DC1280		High value	1.87	± 0.56

【NOTE】

1. The control should be stored in the refrigerator at 2°C ~8°C. After reconstitution, the control can be stable for 8 hours at room temperature (10°C~30°C), and can be stable for 12 hours at 2°C~8°C.
2. The control cannot be frozen after reconstitution.
3. The control is for IVD only.