




# Table of C-Reactive Protein Calibrator (CRP) Reference Value

## C-Reactive Protein Calibrator

LOT 2025073101

 2026-07-30

Applicable Instruments	DH51CRP, DH53CRP, DH56CRP, D1-CRP, D3-CRP, D5-CRP (Technical file version A11.10 or higher with heating)	DH51CRP, DH53CRP, DH56CRP, D1-CRP, D3-CRP, D5-CRP (Technical file version A11.10 or higher without heating)	DH71CRP, DH73CRP, D2-CRP, D7-CRP (Technical file version A6.0 or higher with heating)	DH71CRP, DH73CRP, D2-CRP, D7-CRP (Technical file version A6.0 or higher without heating)	D6-CRP, DH76CRP (Technical file version B2.0)	D6-CRP, DH76CRP (Technical file version B2.2 or higher)	DF1-CRP, DF3-CRP, DF5-CRP, DF50CRP, DF52CRP, DF53CRP (Technical file version A6.0 or higher with heating)
Level (multi-point)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)
Calibrator 1	6.3	6.6	5.3	6.3	5.4	5.3	5.6
Calibrator 2	21.4	20.9	21.6	22.7	19.7	20.2	19.4
Calibrator 3	42.2	42.5	42.6	42.1	43.0	39.5	39.5
Calibrator 4	131.8	110.5	111.9	113.2	103.5	95.2	141.3
Calibrator 5	177.7	155.8	149.3	145.4	128.4	118.5	166.2
Calibrator 6	180.4	175.0	168.3	173.8	158.0	160.2	182.6
Applicable Instruments	DF1-CRP, DF3-CRP, DF5-CRP, DF50CRP, DF52CRP, DF53CRP (Technical file versions A8.1 and below without heating)	DM79X, DM78X, DM72X (Technical file version A1.0 or higher)	DH36X, DH33X, DH31X (Technical file versions A1.0 or higher)	DP-H10, DP-H12, DP-H13, DP-H15, DP-H16, DP-H17 (Technical file versions A1.0 or higher)	DH-800 CRP series DH-800 CS series (Technical file Version A1.3 or higher)	DH-2000T series , DH-100T, DH-102T DH-103T, DH-105T DH-106T, DH-108T (Technical file version A2.0 or higher)	/
Level (multi-point)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	Reference Value (mg/L)	/
Calibrator 1	5.9	5.95	5.4	5.6	4.44	4.91	/
Calibrator 2	24.9	22.42	21.0	24.3	19.55	20.55	/
Calibrator 3	45.9	44.93	42.2	52.7	38.42	40.80	/
Calibrator 4	122.8	100.92	108.8	111.3	80.39	93.17	/
Calibrator 5	155.4	132.15	148.8	153.4	102.48	128.97	/
Calibrator 6	190.7	170.61	172.5	184.6	165.76	162.88	/