




# Table of Serum Amyloid A Protein Control (SAA) Reference Value

## Serum Amyloid A Protein Control

LOT 2025081201

 2026-08-11

Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)	Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)
DM79X、DM78X、DM72X (Applicable to technical file version A1.0 or higher)	Control 1	8.98	±2.69	DP-H10、DP-H12、DP-H13、DP-H15、DP-H16、DP-H17 (Applicable to technical file version A1.0 or higher)	Control 1	8.0	±2.40
	Control 2	52.72	±15.82		Control 2	57.8	±17.34