


Table of Serum Amyloid A Protein Control (SAA) Reference Value

Serum Amyloid A Protein Control

LOT 2025082601

 2026-08-25

| Applicable Instruments | Level | Reference Value (mg/L) | Deviation Range (mg/L) | Applicable Instruments | Level | Reference Value (mg/L) | Deviation Range (mg/L) |
|---------------------------------------------------------------------------------------------------------------------|-----------|------------------------|------------------------|-------------------------------------------------------------|-----------|------------------------|------------------------|
| DH-2000T series DH-100T、DH-102T DH-103T、DH-105T DH-106T、DH-108T (Technical File version A2.0 or higher) | Control 1 | 8.66 | ±2.60 | DH-800 CS series (Technical File Version A1.3 or higher) | Control 1 | 8.03 | ±2.41 |
| | Control 2 | 46.38 | ±13.91 | | Control 2 | 55.49 | ±16.65 |