




|  | 文件类型 | 技术文件 | 文件版本 | V3.6.0 | |
|---|---|-------------------------|------|------------------------|------------|
| | 项目代号 | / | 页码 | 第0页 共1页 | |
| | 文件编号 | GL DM-56-62-en | | | |
| 文件名称： <p style="text-align: center;">CONTROL DM-BF BF1125 2026-01-10</p> 适用范围： 1117项目 | | | | | |
| 编制：刘金盈 日期：2025.10.17 | | 审核：冯会娟 日期：2025.10.17 | | 批准：程军 日期：2025.10.17 | |
| 修 订 记 录 | | | | | |
| 版本 | TCN/ECR/PCN编号 | 修订内容概要 | 修改人 | 批准人 | 批准日期 |
| V3.6.0 | / | 新建 | 刘金盈 | 程军 | 2025.10.17 |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| 发至： | <input type="checkbox"/> 总经办 <input type="checkbox"/> 质量部 <input type="checkbox"/> 财务部 <input checked="" type="checkbox"/> 营销部 <input checked="" type="checkbox"/> 客服 <input type="checkbox"/> 计划 <input type="checkbox"/> 采购 <input checked="" type="checkbox"/> 工程 <input type="checkbox"/> 仓库 <input checked="" type="checkbox"/> 设备研发 <input checked="" type="checkbox"/> 试剂研发 <input type="checkbox"/> 物料质量 <input type="checkbox"/> 设备生产 <input type="checkbox"/> 试剂生产 <input checked="" type="checkbox"/> 试剂质检 <input checked="" type="checkbox"/> 设备质检 | | | | |
| 由工程 | | | | | |

DM-BF

HEMATOLOGY CONTROL

Reference Values provided by DYMIND
CONTROL
 2025-10-13

 2026-01-10

| Applicable Instruments | Parameter | Unit | LOT BF1125L | LOT BF1125N | LOT BF1125H |
|---|-----------|---------------------|---|---|---|
| DYMIND DH-800 series DH-800 CRP series DH-800 CS series (Technical File Version A1.3 or higher) | WBC-BF | 10 ⁹ /L | 0.125 ±0.050 | 0.266 ±0.150 | 0.912 ±0.200 |
| | RBC-BF | 10 ¹² /L | 0.029 ±0.020 | 0.080 ±0.060 | 0.480 ±0.100 |

【NOTE】

1. The controls should be stored in refrigerator (2℃~8℃). After opening, it will keep stable for 30 days when it is stored airtight at 2℃~8℃.
2. Before mixing and running the control after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15℃~30℃).
- 3.Controls must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
4. After using, put the controls back into the refrigerator to prevent contamination and evaporation.



BF1125L



BF1125N



BF1125H