


Fibrin and Fibrinogen Degradation Product Calibrator


Reference Value Table

FDP Calibrator

CALIBRATOR

LOT 2025062901

 20250629

 20270628

Applicable Instruments	Unit	Reference Value
CA1200	μg/mL	FDP calibrator 1: 4.75
CA1203		FDP calibrator 2: 20.64
CA1205		FDP calibrator 3: 40.39
DC1210		FDP calibrator 4: 73.82
DC1250		FDP calibrator 5: 107.75
DC1280		

【NOTE】

1. The calibrator should be stored in the refrigerator at 2 °C ~ 8 °C. After reconstitution, the calibrator can be stable for 8 hours at room temperature (10 °C~ 30 °C), and can be stable for 12 hours at 2 °C ~ 8 °C.
2. The calibrator cannot be frozen after reconstitution.
3. The calibrator is for IVD only.