


Table of C-Reactive Protein Control (CRP) Reference Value

C-Reactive Protein Control

LOT 2025091801


 2026-09-17

Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)	Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)
DH51CRP, DH53CRP, DH56CRP, D1-CRP, D3-CRP, D5-CRP (Technical file version A11.10 or higher with heating)	Control 1	3.3	±1.0	DH51CRP, DH53CRP, DH56CRP, D1-CRP, D3-CRP, D5-CRP (Technical file version A11.10 or higher without heating)	Control 1	3.2	±1.0
	Control 2	14.1	±4.2		Control 2	14.3	±4.3
	Control 3	33.7	±10.1		Control 3	33.6	±10.1
DH71CRP, DH73CRP, D2-CRP, D7-CRP (Technical file version A6.0 or higher with heating)	Control 1	2.7	±0.8	DH71CRP, DH73CRP, D2-CRP, D7-CRP (Technical file version A6.0 or higher without heating)	Control 1	2.8	±0.8
	Control 2	14.9	±4.5		Control 2	13.7	±4.1
	Control 3	39.4	±11.8		Control 3	35.7	±10.7
DH-2000T series , DH-100T, DH-102T, DH-103T, DH-105T, DH-106T, DH-108T (Technical file version A2.0 or higher)	Control 1	2.84	±0.85	D6-CRP, DH76CRP (Technical file version B2.2 or higher)	Control 1	2.5	±0.8
	Control 2	14.53	±4.36		Control 2	12.9	±3.9
	Control 3	36.61	±10.98		Control 3	34.4	±10.3
DF1-CRP, DF3-CRP, DF5-CRP, DF50CRP, DF52CRP, DF53CRP (Technical file version A6.0 or higher with heating)	Control 1	2.5	±0.8	DF1-CRP, DF3-CRP, DF5-CRP, DF50CRP, DF52CRP, DF53CRP (Technical file versions A8.1 and below without heating)	Control 1	3.3	±1.0
	Control 2	11.8	±3.5		Control 2	15.5	±4.6
	Control 3	33.4	±10.0		Control 3	38.2	±11.5

Table of C-Reactive Protein Control (CRP) Reference Value

C-Reactive Protein Control

LOT 2025091801

 2026-09-17

Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)	Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)
DM79X, DM78X, DM72X (Technical file version A1.0 or higher)	Control 1	3.63	±1.1	DH36X, DH33X, DH31X (Technical file versions A1.0 or higher)	Control 1	2.7	±0.8
	Control 2	15.97	±4.8		Control 2	12.5	±3.8
	Control 3	40.04	±12.0		Control 3	34.9	±10.5
DH-800 CRP series DH-800 CS series (Technical file Version A1.3 or higher)	Control 1	3.06	±0.92	DP-H10, DP-H12 DP-H13, DP-H15 DP-H16, DP-H17 (Technical file versions A1.0 or higher)	Control 1	4.1	±1.2
	Control 2	15.69	±4.71		Control 2	19.0	±5.7
	Control 3	38.12	±11.44		Control 3	45.6	±13.7