



Fibrinogen Calibrator Reference Value Table

Fibrinogen Calibrator

CALIBRATOR

LOT 2025101601

 20251016 20271015

Applicable Instrument	Unit	Principle	Reference Value
CA1200 CA1203 CA1205 DC1210	g/L	Optical method	2.81
DC1250 DC1280	g/L	Magnetic bead	2.89

【NOTE】

1. The calibrator should be stored in the refrigerator at 2 °C - 8 °C. After reconstitution, the calibrator can be stable for 4 hours at room temperature (10 °C - 30 °C), and can be stable for 8 hours at 2 °C - 8 °C.
2. The calibrator cannot be frozen after reconstitution
3. The calibrator is for IVD only.