




# Table of Serum Amyloid A Protein Control (SAA) Reference Value

## Serum Amyloid A Protein Control

LOT 2025110101

 2026-10-31

Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)	Applicable Instruments	Level	Reference Value (mg/L)	Deviation Range (mg/L)
DM79X、DM78X、DM72X (Applicable to technical file version A1.0 or higher)	Control 1	9.67	±2.90	DP-H10、DP-H12、DP-H13、DP-H15、DP-H16、DP-H17 (Applicable to technical file version A1.0 or higher)	Control 1	8.0	±2.40
	Control 2	49.13	±14.74		Control 2	54.0	±16.20