


Table of Serum Amyloid A Protein Calibrator (SAA) Reference Value

Serum Amyloid A Protein Calibrator

LOT 2025110101

 2026-10-31

Applicable Instruments	DM79X、DM78X、DM72X (Applicable to technical file version A1.0 or higher)	DP-H10、DP-H12、DP-H13、DP-H15、DP-H16、 DP-H17 (Applicable to technical file version A1.0 or higher)
Level (multi-point)	Reference Value (mg/L)	Reference Value (mg/L)
Calibrator 1	11.00	9.6
Calibrator 2	19.20	19.2
Calibrator 3	49.13	54.0
Calibrator 4	95.00	100.3
Calibrator 5	169.00	191.0
Calibrator 6	189.00	/
Calibrator 7	229.00	/